

Jazzy 1143

MIGHTY


The Jazzy® 1143 Mighty is designed for excellent indoor maneuverability and outstanding outdoor capability. With 12" drive wheels for enhanced outdoor performance and a 300 lbs. weight capacity, combined with patented mid-wheel drive performance and maneuverability, the Jazzy 1143 Mighty is the ultimate indoor/outdoor power chair.



Jazzy 1143
shown in Onyx Black
with gray contoured
medium-back seat,
foot platform and
VSI controller

Weight capacity:	300 lbs.
Maximum speed:	Up to 4.5 mph
Ground clearance:	4"
Turning radius:	20.5"
Overall length:	39"
Overall width:	25"
Wheels:	
Drive:	12" pneumatic
Casters:	9" pneumatic
Anti-tips:	8" solid
Suspension:	Limited
Drive train:	Two-motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Electronics:	50 amp, PG VSI controller
Specialty controls:	Not accepted
Battery charger:	On-board, 4 amp (standard)
Per-charge range:	Up to 25 miles
Battery requirements:	(2) 12 volt deep cycle
Battery size:	U-1, 32 AH
Transport weight:	
Battery weight:	24.5 lbs. each
Base weight:	97 lbs.
Std. seat weight:	33 lbs. (medium-back)
Seat dimensions:	20" x 20" (max. contoured seating)
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors

Colors:



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

FEATURES

- Powerful two-motor, mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 12" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4.5 mph
- On-board battery charger
- Weight capacity of 300 lbs.
- Numerous seating options including the Synergy™ seating system
- 8" resilient front anti-tips
- Fixed-height, flip-back armrests

OPTIONS

- Contoured high-back seat
- Synergy™ seating system
- Power elevating seat
- Contoured medium/high-back with solid seat pan
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Flat-free tires
- Adjustable-height armrests (optional on standard seats)
- Removable height-adjustable arms (optional on Synergy Seat)
- Swing-away joystick mount
- Safety flag