

QUANTUM[®]

1121



The Quantum 1121 with Active-Trac[®] Suspension features powerful new GLIDE motors for added power & torque plus a narrow overall width of only 23.8" making it an extremely maneuverable, high performance power chair solution indoors and out.




1121

SPECIFICATIONS

DRIVE WHEELS:	14" pneumatic
REAR CASTERS:	8" solid
ANTI-TIPS:	6" solid, adjustable, front-mounted
SUSPENSION:	Active-Trac®
MAXIMUM SPEED:	Up to 5 mph ¹
BRAKES:	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE:	3.5"
TURNING RADIUS:	20.5"
OVERALL SIZE:	Length: 38.5" Width: 23.8"
SEATING:	Synergy® Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® Power Positioning System Adjustable Width: 14"-20" Adjustable Depth: 14"-20" Specialty Seat (standard) Power elevating seat
DRIVE TRAIN:	Two-motor, mid-wheel drive
BATTERIES:	Two 12 volt deep cycle NF-22
RANGE PER CHARGE:	Up to 25 miles ¹
BATTERY CHARGER:	Lester 8A, Off-board
ELECTRONICS:	70A PG VSI
WEIGHT CAPACITY:	300 lbs. maximum (user weight)
BASE WEIGHT:	139.5 lbs.
SEAT WEIGHT:	37 lbs. Specialty Seat
BATTERY WEIGHT:	38 lbs. (each)
COLOR CHOICES:	Candy Apple Red, Viper Blue, Onyx Black, Champagne, Sunburst Orange
WARRANTY:	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor

¹) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.

AVAILABLE COLORS



FEATURES

- Powerful new GLIDE motors for added power & torque
- Active-Trac® Suspension for enhanced performance on rough and uneven terrain, creating a smooth and stable ride
- Rear articulating caster beam for great performance and durability
- One-piece, 2" diameter tubular frame for maximum strength and enhanced performance
- Easy front access to batteries
- 70A PG VSI controller
- Standard Specialty Seat with swing-away leg rests

OPTIONS

- Contoured seating
- Solid tires
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Rear basket
- Angle-adjustable footplates
- 70° and 80° heavy-duty drop-in footrests
- On-board/Off-board charger option

QUANTUM 1121

Shown in Viper Blue with Synergy® Seating System and Remote Plus controller



QUANTUM
REHAB
INNOVATIVE REHAB SOLUTIONS
A Division of Pride Mobility Products® Corp.

1-866-800-2002 (US) • 1-888-570-1113 (Canada) • quantumrehab.com

It is recommended that you consult a qualified Rehab Provider prior to fitting or ordering any Quantum Rehab product. The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.