

QUANTUM[®]

1402



Quantum[®] 1402
Shown in Candy Apple Red
with Synergy[®] Seat and
PG VSI controller

The Quantum[®] 1402 delivers superb rehab capability in a compact, sporty design that is available with a choice of motor packages for unmatched versatility: a standard high-torque package with a 300 lbs. weight capacity and a speed of 4.5 mph; an optional 6 mph high-speed package with 300 lbs. weight capacity; or an optional Hammer motor package with a weight capacity of 400 lbs. and increased speed of up to 6.5 mph.



www.quantumrehab.com

Q1402

SPECIFICATIONS	
DRIVE WHEELS:	14" pneumatic
REAR CASTERS:	8" solid
ANTI-TIPS:	6" solid, adjustable, front-mounted
SUSPENSION:	Active-Trac [®] with rear caster suspension
STANDARD SPEED:	Up to 4.5 mph ¹ (high-torque motors) ⁵
BRAKES:	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE:	3.25"
TURNING RADIUS:	21"
OVERALL SIZE:	Length: 40.5" Width: 25.4"
SEATING:	Synergy™ Seat: Adjustable Width: 10"-24" Adjustable Depth: 10"-24" TRU-Balance™ Super-Low Power Tilt ² TRU-Balance™ Power Tilt (available with manual recline): Adjustable Width: 14"-24" Adjustable Depth: 14"-24" Contoured high-back Contoured medium-back Manual recliner Solid seat pan option Power elevating seat
BATTERIES:	Two 12 volt, deep cycle Group 24
RANGE PER CHARGE:	Up to 25 miles ¹
BATTERY CHARGER²:	8A, Off-board
ELECTRONIC:	70A, PG VSI controller (standard) ³
WEIGHT CAPACITY:	300 lbs. maximum (user weight) ⁴
BASE WEIGHT:	140 lbs.
SEAT WEIGHT:	37 lbs. (medium-back)
BATTERY WEIGHT:	53.5 lbs. (each)
COLOR CHOICES:	Candy Apple Red, Viper Blue, Onyx Black and Champagne, Sunburst Orange
WARRANTY:	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor

¹) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.

²) TRU-Balance™ Tilt / super low tilt requires off-board charger.

³) Hammer motor package requires 100 amp Remote Plus controller.

⁴) 400 lbs. weight capacity when Hammer motor package selected.

⁵) See Options box for other speed configurations.

FEATURES
<ul style="list-style-type: none"> • Easy front-access batteries • Powerful, patented two-motor mid-wheel drive design provides superb power and acceleration • Active-Trac[®] and rear caster suspension for outstanding outdoor performance • Adjustable front anti-tips to meet your preferences and terrain challenges • Battery-conserving sleep mode available with Remote Plus controller

OPTIONS
<ul style="list-style-type: none"> • High-speed motor package (top speed of 6 mph with a weight capacity of 300 lbs.) • Hammer motor package (top speed of 6.5 mph with a weight capacity of 400 lbs.)³ • TRU-Balance™ power positioning seating system • Residual limb supports • Solid tires • Lap belt • Cup holder • Heavy-duty dust cover • Swing-away joystick • Swing-away footrests • Elevating leg rests • Oxygen tank holder • Cane/crutch holder • Multi-axis joystick • Saddlebag • Walker holder • Rear basket • Specialty controls • Angle-adjustable foot plates • 70° and 80° heavy-duty, drop-in footrests • Adductor buttons • Neoprene footrest sleeves • On board/Offboard charger combo

AVAILABLE COLORS
    
<p>Candy Apple Red Viper Blue Onyx Black Champagne Sunburst Orange</p>