

QUANTUM[®]

1420



QUANTUM 1420

Shown in Viper Blue with Synergy™ Seat and PG Remote Plus controller

The Quantum 1420 delivers remarkable rehab versatility along with exclusive performance features like Active-Trac[®] Suspension, two-motor mid-wheel drive design, optional high-speed motor package, and adjustable front anti-tips for enhanced outdoor capability and superb tight-quarter maneuverability. This combination of performance and versatility makes the Quantum 1420 the ideal power chair for all your rehab needs.



www.quantumrehab.com


Q1420

SPECIFICATIONS	
DRIVE WHEELS:	16" pneumatic
REAR CASTERS:	9" pneumatic
ANTI-TIPS:	8" solid, adjustable, front-mounted
SUSPENSION:	Active-Trac® with rear caster suspension
MAXIMUM SPEED:	Up to 4.5 mph ¹
BRAKES:	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE:	4.5"
TURNING RADIUS:	22.5"
OVERALL SIZE:	Length: 45.5" Width: 26.75"
SEATING:	Synergy™ Seat: Adjustable Width: 10"-24" Adjustable Depth: 10"-24" TRU-Balance™ Super-Low Power Tilt ² TRU-Balance™ Power Tilt (available with manual recline): Adjustable Width: 14"-24" Adjustable Depth: 14"-24" Contoured high-back Contoured medium-back Manual recliner Solid seat pan option Power elevating seat
DRIVE TRAIN:	Two-motor, mid-wheel drive; high speed 6mph (optional) ²
BATTERIES:	Two 12 volt, deep cycle 70 AH, Group 24
RANGE PER CHARGE:	Up to 25 miles ¹
BATTERY CHARGER:	8A, Off-board
ELECTRONICS:	70A, Remote Plus Controller
WEIGHT CAPACITY:	500lbs. maximum (user weight)
BASE WEIGHT:	143 lbs.
SEAT WEIGHT:	41 lbs. (Synergy)
BATTERY WEIGHT:	53.5 lbs. (each)
COLOR CHOICES:	Candy Apple Red, Viper Blue, Onyx Black, Champagne and Sunburst Orange
WARRANTY:	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor

¹) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
²) 425 lbs. weight capacity

FEATURES
<ul style="list-style-type: none"> • High-torque motor package delivers a top speed of 4.5 mph and a weight capacity of 500lbs. • Powerful, patented two-motor mid-wheel drive design provides superb power and acceleration • Active-Trac® and rear caster suspension for outstanding outdoor performance • Adjustable front anti-tips to meet your preferences and terrain challenges • 8A powerful off-board charger • Battery-conserving sleep mode available with Remote Plus controller

OPTIONS
<ul style="list-style-type: none"> • Synergy™ manual recline • Synergy™ manual tilt • TRU-Balance™ Power Positioning Systems • Solid tires • Lap belt • Cup holder • Heavy duty dust cover • Swing-away joystick • Swing-away footrests • Elevating leg rests • Oxygen tank holder • Cane/crutch holder • Walker holder • Rear basket • Specialty controls • Angle-adjustable footplates • 70° and 80° heavy duty drop-in footrests • Saddlebag • High-speed motor package-top speed of 6 mph / 425 lbs. weight capacity • On-board/ Off-board charger combo

AVAILABLE COLORS

<p>Candy Apple Red Viper Blue Onyx Black Champagne Sunburst Orange</p>



1-866-800-2002 (US) • 1-888-570-1113 (Canada) • www.quantumrehab.com

It is recommended that you consult a qualified Rehab Provider prior to fitting or ordering any Quantum Rehab product. The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.