

QUANTUM®

600

&

QUANTUM®

600 XL

(Quantum 600 Upgrade)



Q6 TORQUE
FEEL THE POWER

QUANTUM® 600
(Available with XL option)
shown in Sunburst Orange
with optional TRU-Balance®
Positioning System and
PG Remote Plus controller

Equipped with our patent pending Mid-Wheel 6™, the Quantum® 600 and Quantum® 600 XL deliver maximum stability without compromising performance. The independent active OMNI-Casters (nylon, spherical-shaped casters) raise and lower to compensate instantly for changes in the terrain and allow for smooth and stable transitions. In addition, Quantum Rehab has equipped the units with powerful new GLIDE motors for added power and torque.



Q6 | POWER CHAIR SERIES

SPECIFICATIONS

DRIVE WHEELS:	14" pneumatic / XL: 14" solid
OMNI-CASTERS:	Front: 5" solid Rear: 6" solid
SUSPENSION:	ATX Suspension (Active-Trac® with extra stability)
MAXIMUM SPEED:	Up to 5 mph ¹ / XL: 4.5 mph
BRAKES:	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE:	3"
TURNING RADIUS:	20" / XL: 21"
OVERALL SIZE:	Length: 37.5" Width: 23.5"
SEATING:	Synergy® Seat Adjustable Width: 10"-20" / XL: 10"-24" Adjustable Depth: 10"-20" / XL: 10"-24" TRU-Balance® Power Positioning System ² Adjustable Width: 14"-20" / XL: 10"-24" Adjustable Depth: 14"-20" / XL: 10"-24"
DRIVE TRAIN:	Two-motor, mid-wheel drive
BATTERIES:	Two 12 volt, deep cycle NF-22
RANGE PER CHARGE:	Up to 25 miles ¹ / XL: Up to 15 miles ¹
BATTERY CHARGER:	8A Lester, off-board (standard)
ELECTRONICS:	75A, Flight controller
WEIGHT CAPACITY:	300 lbs. maximum (user weight) / XL: 400 lbs.
BASE WEIGHT:	130 lbs. / XL: 136.5 lbs.
SEAT WEIGHT:	37 lbs. (medium-back)
BATTERY WEIGHT:	38 lbs. (each)
WARRANTY:	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor

COLORS:



Candy
Apple Red



Viper Blue



Onyx
Black



Champagne



Sunburst
Orange

- ¹) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- ²) TRU-Balance® Power Recline not compatible with XL model.

FEATURES

- Mid-Wheel 6™ allows six wheels on the ground for maximum stability
- OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters for enhanced performance over more varied terrain
- 4-pole Hammer motors

OPTIONS

- Swing-away joystick
- Multi-axis joystick mount
- Elevating leg rests
- Oxygen tank holder
- Neoprene leg rest sleeves
- Gel adductor buttons
- Angle-adjustable footplates
- 70° drop-in footrests
- Lap belt
- On/off-board charger combo
- High-back seat



888-570-1113 • St. Catharines, Ontario
Visit us on the web at www.quantumrehab.com