

# QUANTUM® 6000

&

# QUANTUM® 6000 XL



**Q6 TORQUE**  
**FEEL THE POWER**

QUANTUM® 6000  
(Available with XL option)  
shown in Metallic Red with TRU-Balance®  
Power Tilt & Recline and Articulating  
Leg Rests

The Quantum® 6000 & Quantum® 6000 XL incorporate aggressive new features, powerful 4-pole Hammer motors, and our patent pending Mid-Wheel 6™ design platform to deliver unmatched performance and tight-space maneuverability while delivering rehab capability.



Q6 | POWER CHAIR SERIES

## SPECIFICATIONS

<b>DRIVE WHEELS:</b>	14" pneumatic / <b>XL:</b> 14" solid
<b>OMNI-CASTERS:</b>	<b>Front:</b> 5" solid <b>Rear:</b> 6" solid
<b>SUSPENSION:</b>	Active-Trac® ATX Suspension with integrated rear dual-strut suspension and articulating caster beam
<b>MAXIMUM SPEED:</b>	Up to 6.5 mph <sup>1</sup>
<b>BRAKES:</b>	Intelligent Braking (electronic, regenerative disc brakes)
<b>GROUND CLEARANCE:</b>	3.5"
<b>TURNING RADIUS:</b>	22"
<b>OVERALL SIZE:</b>	Length: 38" (without front riggings) Width: 25.75"
<b>OPTIONAL SEATING:</b>	Synergy® Seat Adjustable Width: 10"-20"/ <b>XL:</b> 10"-24" Adjustable Depth: 10"-20"/ <b>XL:</b> 10"-24" TRU-Balance® Power Positioning System <sup>2</sup> Adjustable Width: 14"-20"/ <b>XL:</b> 14"-24" Adjustable Depth: 14"-20"/ <b>XL:</b> 14"-24" High-back (max. dimensions: 22"x20") (max dimensions <b>XL:</b> 24"x24") Medium-back (max. dimensions: 22"x20") (max. dimensions <b>XL:</b> 24"x24")
<b>DRIVE TRAIN:</b>	Two-motor, mid-wheel drive
<b>BATTERIES:</b>	Two 12 volt, deep cycle Group 24
<b>RANGE PER CHARGE:</b>	Up to 25 miles <sup>1</sup> / <b>XL:</b> Up to 15 miles <sup>1</sup>
<b>BATTERY CHARGER:</b>	8A, Lester off-board
<b>ELECTRONICS:</b>	100A, PG Remote Plus controller
<b>WEIGHT CAPACITY:</b>	300 lbs. maximum (user weight)/ <b>XL:</b> 400 lbs.
<b>BASE WEIGHT:</b>	168.5 lbs.
<b>SEAT WEIGHT:</b>	37 lbs. (Specialty Seat)
<b>BATTERY WEIGHT:</b>	53.5 lbs. (each)
<b>WARRANTY:</b>	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor

### COLORS:



Metallic Red

Metallic Blue

Metallic Silver

Sunburst Orange

Metallic Black

<sup>1</sup>) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.

<sup>2</sup>) TRU-Balance® Power Recline not compatible with XL model.

## FEATURES

- Mid-Wheel 6™ allows six wheels on the ground for maximum stability
- OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters for enhanced performance over more varied terrain
- 4-pole Hammer motors
- Accepts (2) Group 24 & (1) NF-22 when used with a vent tray

## OPTIONS

- Swing-away joystick
- Multi-axis joystick mount
- Elevating leg rests
- Oxygen tank holder
- Neoprene leg rest sleeves
- Gel adductor buttons
- Vent tray
- Angle-adjustable footplates
- Lap belt
- Contoured seating
- High-mount foot platform
- 70° HD drop-in footrests
- Residual limb supports
- Contoured seating with solid seat pan
- Power elevating seat



888-570-1113 • St. Catharines, Ontario  
Visit us on the web at [www.quantumrehab.com](http://www.quantumrehab.com)