

# QUANTUM<sup>®</sup> 600E

**Q6 TORQUE**  
**FEEL THE POWER**



**QUANTUM<sup>®</sup> 600E**  
shown in Candy Apple Red with Specialty  
Seat Kit and PG VSI controller

The Quantum<sup>®</sup> 600E is the perfect balance of performance and value. The independent active OMNI-Casters (nylon, spherical-shaped casters) raise and lower to compensate instantly for changes in terrain and allow for smooth and stable transitions. In addition, Quantum<sup>®</sup> Rehab has equipped the unit with new powerful motors for added power and torque, and a standard off-board charger is featured.

**Q6**  
Q6 | POWER CHAIR SERIES

## SPECIFICATIONS

<b>DRIVE WHEELS:</b>	14" pneumatic
<b>OMNI-CASTERS:</b>	<b>Front:</b> 5" solid <b>Rear:</b> 6" solid
<b>SUSPENSION:</b>	ATX Suspension (Active-Trac® with extra stability)
<b>MAXIMUM SPEED:</b>	Up to 5 mph <sup>1</sup>
<b>BRAKES:</b>	Intelligent Braking (electronic, regenerative disc brakes)
<b>GROUND CLEARANCE:</b>	3"
<b>TURNING RADIUS:</b>	20"
<b>OVERALL SIZE:</b>	Length: 36.5" Width: 23.5"
<b>SEATING:</b>	Synergy® Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® Power Positioning System Adjustable Width: 14"-20" Adjustable Depth: 14"-20" Medium-back (standard)
<b>DRIVE TRAIN:</b>	Two-motor, mid-wheel drive
<b>BATTERIES:</b>	Two 12 volt, deep cycle NF-22
<b>RANGE PER CHARGE:</b>	Up to 25 miles <sup>1</sup>
<b>BATTERY CHARGER:</b>	8A Lester, off-board (standard)
<b>ELECTRONICS:</b>	70A, PG VSI controller
<b>WEIGHT CAPACITY:</b>	300 lbs. maximum (user weight)
<b>BASE WEIGHT:</b>	130 lbs.
<b>SEAT WEIGHT:</b>	37 lbs. (medium-back)
<b>BATTERY WEIGHT:</b>	38 lbs. (each)
<b>WARRANTY:</b>	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor

### COLORS:



Candy  
Apple Red



Viper Blue



Onyx  
Black



Champagne



Sunburst  
Orange

<sup>1</sup>) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.

## FEATURES

- Mid-Wheel 6™ allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® Power Positioning
- Patent pending OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- Patent pending ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters for enhanced performance over more varied terrain
- Easy front access to batteries
- 70A, PG VSI controller
- Powerful new motors for added power and torque
- Off-board charger

## OPTIONS

- Solid tires
- Swing-away joystick
- Multi-axis joystick mount
- Elevating leg rests
- Oxygen tank holder
- Neoprene leg rest sleeves
- Gel adductor buttons
- Angle-adjustable footplates
- 70° drop-in footrests
- Lap belt
- High-mount foot platform
- Synergy® Seat
- Cane/crutch holder
- Specialty Seat

**QUANTUM**  
**REHAB**  
INNOVATIVE REHAB SOLUTIONS  
A Division of Pride Mobility Products® Corp.

866-800-2002 (US) ▪ 888-570-1113 (Canada)

Visit us on the web at [www.quantumrehab.com](http://www.quantumrehab.com)