



Q6 | POWER CHAIR SERIES

# QUANTUM<sup>®</sup> 610



FEATURING  
MID • WHEEL



TM

**QUANTUM<sup>®</sup> 610**  
shown in Candy Apple Red  
with optional Synergy<sup>®</sup> Seating  
System and PG Remote Plus  
controller

The Quantum<sup>®</sup> 610 combines a compact frame with our patent pending Mid-Wheel 6<sup>TM</sup> design platform to provide the ultimate choice for users who demand maximum stability, tight maneuverability and true all-around rehab performance. In addition, Quantum Rehab has equipped the unit with powerful new GLIDE motors for added power and torque.

Find out more about the Quantum<sup>®</sup> 610 on the web at [www.quantumrehab.com](http://www.quantumrehab.com)

## SPECIFICATIONS

<b>DRIVE WHEELS:</b>	10" pneumatic
<b>OMNI-CASTERS:</b>	<b>Front:</b> 5" solid <b>Rear:</b> 6" solid
<b>SUSPENSION:</b>	ATX Suspension (Active-Trac® with extra stability)
<b>MAXIMUM SPEED:</b>	Up to 4 mph <sup>1</sup>
<b>BRAKES:</b>	Intelligent Braking (electronic, regenerative disc brakes)
<b>GROUND CLEARANCE:</b>	2.125"
<b>TURNING RADIUS:</b>	21.25"
<b>OVERALL SIZE:</b>	Length: 35.5" Width: 23.5"
<b>SEATING:</b>	Synergy® Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® Power Positioning System <sup>2</sup> Adjustable Width: 14"-20" Adjustable Depth: 14"-20" Specialty Seat Kit High-back Medium-back
<b>DRIVE TRAIN:</b>	Two-motor, mid-wheel drive
<b>BATTERIES:</b>	Two 12 volt, deep cycle U-1
<b>RANGE PER CHARGE:</b>	Up to 25 miles <sup>1</sup>
<b>BATTERY CHARGER:</b>	5A Lester, PowerCube off-board (standard)
<b>ELECTRONICS:</b>	75A Flight Controller
<b>WEIGHT CAPACITY:</b>	300 lbs. maximum (user weight)
<b>BASE WEIGHT:</b>	126.5 lbs.
<b>BATTERY WEIGHT:</b>	24.5 lbs. (each)
<b>WARRANTY:</b>	Lifetime limited warranty on frame 2-year warranty on electronics 18-month warranty on drive motor

### COLORS:



Candy  
Apple Red



Viper Blue



Onyx  
Black



Champagne



Sunburst  
Orange

## FEATURES

- Mid-Wheel 6™ allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® Power Positioning
- OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters for enhanced performance over more varied terrain
- Easy front access to batteries
- 75A Flight Controller
- Powerful new GLIDE motors for added power and torque
- Off-board charger (5A PowerCube)

## OPTIONS

- Solid tires
- Swing-away joystick
- Multi-axis joystick mount
- Elevating leg rests
- Oxygen tank holder
- Neoprene leg rest sleeves
- Gel adductor buttons
- Power elevating seat
- Angle-adjustable footplates
- 70° HD drop-in footrests
- Lap belt
- On/off-board charger combo

### FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) 200 lbs. weight capacity with TRU-Balance® Power Positioning System; off-board charger required

# QUANTUM REHAB

INNOVATIVE REHAB SOLUTIONS

A Division of Pride Mobility Products® Corp.

888-570-1113 • St. Catharines, Ontario • [www.quantumrehab.com](http://www.quantumrehab.com)

It is recommended that you consult a qualified Rehab Provider prior to fitting or ordering any Quantum Rehab product. The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.