



1103 ULTRA



1107



1121



1122

<b>DRIVE WHEELS:</b>	10" pneumatic <sup>1</sup>	10" pneumatic <sup>1</sup>	14" pneumatic <sup>1</sup>	14" pneumatic <sup>1</sup>
<b>CASTERS:</b>	8" solid	6" solid	14" solid	14" solid
<b>ANTI-TIPS:</b>	5" solid	5" solid	8" solid	8" solid
<b>SUSPENSION:</b>	Active-Trac® Suspension	Advanced Active-Trac® Suspension <sup>2</sup>	Active-Trac® Suspension	Active-Trac® Suspension <sup>2</sup>
<b>MAXIMUM SPEED:</b>	Up to 4 mph <sup>6</sup>	Up to 4.5 mph <sup>6</sup>	Up to 5 mph <sup>6</sup>	Up to 6 mph <sup>3,6</sup>
<b>BRAKES:</b>	Electronic, regenerative disc brakes	Electronic, regenerative disc brakes	Electronic, regenerative disc brakes	Electronic, regenerative disc brakes
<b>GROUND CLEARANCE:</b>	2.5"	2.5"	3.5"	3.25"
<b>TURNING RADIUS:</b>	18.5"	20.5"	20.5"	21"
<b>OVERALL SIZE:</b>	Length: 35.5" Width: 23"	Length: 36.5" Width: 23"	Length: 38.5" Width: 23.8"	Length: 40.5" Width: 25.4"
<b>DRIVE TRAIN:</b>	Two-motor, mid-wheel drive	Two-motor, mid-wheel drive	Two-motor, mid-wheel drive	Two-motor, mid-wheel drive
<b>BATTERIES:</b>	Two 12 volt, U-1	Two 12 volt, U-1	Two 12 volt, NF-22	Two 12 volt, Group 24
<b>RANGE PER CHARGE:</b>	Up to 25 miles <sup>6</sup>	Up to 25 miles <sup>6</sup>	Up to 25 miles <sup>6</sup>	Up to 25 miles <sup>6</sup>
<b>BATTERY CHARGER:</b>	4A, On-board <sup>7</sup>	5A, Lester Off-board	8A, Lester Off-board	8A, Lester Off-board
<b>ELECTRONICS:</b>	50A, 6-key VSI <sup>4</sup>	75A, Pride® Flight <sup>4</sup>	70A, PG VSI <sup>4</sup>	70A, PG VSI <sup>4</sup>
<b>WEIGHT CAPACITY:</b>	300 lbs. <sup>5</sup>	300 lbs. <sup>5</sup>	300 lbs. <sup>5</sup>	300 lbs.
<b>BASE WEIGHT:</b>	107 lbs.	112 lbs.	139.5 lbs.	140 lbs.
<b>BATTERY WEIGHT:</b>	24.5 lbs. (each)	24.5 lbs. (each)	38 lbs. (each)	53.5 lbs. (each)

1) Solid tires optional

2) With rear caster suspension

3) With high-speed motor package

4) Available with 70A Remote Plus electronics

5) 200 lbs. weight capacity with TRU-Balance® Power Positioning Systems

6) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.

7) Requires off-board charger with TRU-Balance® Power Positioning System

**COLOR CHOICES:**



Candy Apple Red



Viper Blue



Champagne



Onyx Black



Sunburst Orange  
(Only on 1107 & 1121)

**QUANTUM**  
**REHAB**  
INNOVATIVE REHAB SOLUTIONS  
A Division of Pride Mobility Products® Corp.

866-800-2002 (US) • 888-570-1113 (Canada)

Visit us on the web at [www.quantumrehab.com](http://www.quantumrehab.com)



QUANTUM® 1103 ULTRA

QUANTUM® 1121

QUANTUM® 1107

QUANTUM® 1122

ADVANCED FUNCTIONALITY • ADVANCED DESIGN

**QUANTUM**®

MID-WHEEL DRIVE POWER CHAIRS

## QUANTUM® 1103 ULTRA

THE QUANTUM® 1103 ULTRA brings sporty style and high-performance features, like a standard power elevating seat and Active-Trac® Suspension, together into one perfectly balanced power chair. Yet, even fully loaded, the 1103 Ultra remains highly compact and maneuverable.

### FEATURES

- 300 lbs. weight capacity
- Power elevating seat provides greater access to objects out of reach
- Active-Trac® Suspension (standard) for outstanding outdoor performance
- Powerful two-motor, patented mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 10" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4 mph
- On-board battery charger
- Lap belt



**QUANTUM® 1103 ULTRA**  
with power elevating seat

## QUANTUM® 1107

THE QUANTUM® 1107 easily disassembles to travel with you wherever you go. Yet, complete with advanced Active-Trac® Suspension and powerful new GLIDE motors for added power and torque, the Quantum 1107 offers the stability and performance of a full-size, rigid frame chair.

### FEATURES

- GLIDE motors for added power and torque
- Disassembles for ease of transportation
- Aluminum frame rails and center locking mechanism provide a very rigid platform for ultimate stability
- Innovative battery boxes and front module eliminate battery and motor harnesses
- Advanced Active-Trac® and rear caster suspension for enhanced performance on rough and uneven terrain, creating a smooth and stable ride
- 75 amp, Flight controller
- Standard 5 amp, Lester Power Cube off-board charger



**QUANTUM® 1107**  
with TRU-Balance®  
Power Tilt

## QUANTUM® 1121

THE QUANTUM® 1121 combines outstanding performance with a sleek design. The 1121 utilizes Pride's patented mid-wheel drive technology and an overall length of only 38.5" for great tight-space indoor maneuverability and a turning radius of only 20.5". Equipped with an articulating rear caster beam, patented resilient front anti-tip wheels, Active-Trac® Suspension, and a top speed of 5 mph, the 1121 is ready to conquer all of your daily challenges.

### FEATURES

- Narrow width of 23.8" with a turning radius of 20.5"
- Patented mid-wheel drive design offers exceptional handling
- Active-Trac® Suspension for outstanding outdoor performance
- Accepts a full range of rehab seating including power tilt and recline
- PG VSI, 70 amp controller offers superb control
- Easy front access to batteries without removing the seat or shroud
- Available with optional upgraded electronics for interface with specialty drive controls



**QUANTUM® 1121**  
with Synergy®  
Seating System

## QUANTUM® 1122

THE QUANTUM® 1122, with sleek styling unlike any other power chair, makes a statement all its own. With its fully adaptable rehab capabilities, including an extensive range of seating and electronics options, the Quantum 1122 allows you to make a personal statement as well.

### FEATURES

- Patented mid-wheel drive design and 14" knobby pneumatic tires on mag wheels provide excellent handling and traction
- Standard with high-speed motors for a top speed of 6 mph (4.5 mph with high-torque motors) — choose the performance that best meets your needs
- Active-Trac® and rear caster suspension for outstanding outdoor performance
- Adjustable front anti-tips to meet your preferences and terrain challenges
- Accepts a full range of rehab seating including power tilt and recline and elevating seat combinations
- Available with optional upgraded electronics for interface with specialty drive controls
- Easy front access to batteries without removing the seat or shroud



**QUANTUM® 1122**  
with TRU-Balance®  
Power Positioning  
System