

Q6000

Q600

Q610

DRIVE WHEELS:	14" pneumatic	14" pneumatic	10" pneumatic
OMNI-CASTERS:	Front: 5" solid Rear: 6" solid	Front: 5" solid Rear: 6" solid	Front: 5" solid Rear: 6" solid
SUSPENSION:	Active-Trac® ATX Suspension with integrated rear dual-strut suspension and articulating caster beam	ATX Suspension (Active-Trac® with extra stability)	ATX Suspension (Active-Trac® with extra stability)
MAXIMUM SPEED:	Up to 6.5 mph ¹	Up to 5 mph ¹	Up to 4 mph ¹
BRAKES:	Intelligent Braking (electronic, regenerative disc brakes)	Intelligent Braking (electronic, regenerative disc brakes)	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE:	3"	3"	2.5"
TURNING RADIUS:	22"	20"	21.25"
OVERALL SIZE:	Length: 38" ³ Width: 25.75"	Length: 36.5" ³ Width: 23.5"	Length: 34.75" ³ Width: 23.25"
SEATING:	Synergy® Seat Adjustable Width: 10"-22" Adjustable Depth: 10"-20" TRU-Balance® Power Positioning Adjustable Width: 14"-20" Adjustable Depth: 14"-20" Specialty Seat Kit High-back (22" x 20" max.) Medium-back (22" x 20" max.)	Synergy® Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® Power Positioning Adjustable Width: 14"-20" Adjustable Depth: 14"-20" Specialty Seat Kit High-back (22" x 20" max.) Medium-back (22" x 20" max.)	Synergy® Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® Power Positioning ² Adjustable Width: 14"-20" Adjustable Depth: 14"-20" Specialty Seat Kit High-back (20" x 20" max.) Medium-back (20" x 20" max.)
DRIVE TRAIN:	Two-motor, Mid-Wheel 6®	Two-motor, Mid-Wheel 6®	Two-motor, Mid-Wheel 6®
BATTERIES:	Two 12 volt, deep cycle Group 24	Two 12 volt, deep cycle NF-22	Two 12 volt, deep cycle U-1
RANGE PER CHARGE:	Up to 25 miles ¹	Up to 25 miles ¹	Up to 25 miles ¹
BATTERY CHARGER:	8A Lester, off-board	5A, off-board	5A Lester PowerCube, off-board
ELECTRONICS:	100A, PG Remote Plus controller	75A, Flight controller	75A, Flight controller
WEIGHT CAPACITY:	300 lbs. maximum (user weight)	300 lbs. maximum (user weight)	300 lbs. maximum (user weight)
BASE WEIGHT:	168.5 lbs.	130 lbs.	112 lbs.
BATTERY WEIGHT:	53.5 lbs. (each)	38 lbs. (each)	24.5 lbs. (each)
COLORS:	Metallic Red, Metallic Blue, Metallic Silver, Metallic Black, Sunburst Orange	Candy Apple Red, Viper Blue, Onyx Black, Champagne, Sunburst Orange	Candy Apple Red, Viper Blue, Onyx Black, Champagne, Sunburst Orange



WARRANTY: Lifetime limited warranty on frame; 2-year warranty on electronics; 18-month warranty on drive motor

- ¹⁾ Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- ²⁾ Off-board charger required
- ³⁾ Without front riggings



866-800-2002 (US) • 888-570-1113 (Canada)
Visit us on the web at www.quantumrehab.com

Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %



Conquer Rough Terrain.

Transition Obstacles.

Avoid Getting Stuck.

Q6 TORQUE

FEEL THE POWER™

The Q6 Power Chair Series is designed to provide robust power when it's needed most. The Q6000 Power Chair, with 4-pole Hammer motors, provides a level of torque and efficiency that surpasses other competitive motors. The Q6 Series delivers superior torque and power that achieves "real world" active user satisfaction.

PERFORMANCE

- Patent pending Mid-Wheel 6® drive allows six wheels on the ground for maximum stability
- Superior climbing capability
- Outstanding reverse transitioning
- Powerful motors for added torque and power

HANDLING

- Highly maneuverable and capable of handling aggressive terrain
- Patent pending Active-Trac® (ATX) Suspension for smooth curb transitions and obstacle climbing

SERVICEABILITY

- Standard off-board chargers
- Front access to batteries without seat removal

REHAB CAPABILITY

- Compatible with all TRU-Balance® Power Positioning Systems



Q6 | POWER CHAIR SERIES



QUANTUM® 6000

The Quantum® 6000 incorporates aggressive new features, powerful 4-pole Hammer motors, and our patent pending Mid-Wheel 6® design platform to deliver unmatched performance and tight-space maneuverability while delivering limitless rehab capability. The Quantum 6000 is compatible with all TRU-Balance® Power Positioning Systems and offers an industry-leading vent configuration that allows for battery storage in the base.



ALSO AVAILABLE: *
Quantum 6000 XL
(400 lbs.
weight capacity)



QUANTUM® 600

Equipped with our patent pending Mid-Wheel 6®, the Quantum® 600 delivers maximum stability without compromising performance. The independent active OMNI-Casters (nylon, spherical-shaped casters) raise and lower to compensate instantly for changes in the terrain and allow for smooth and stable transitions. In addition, Quantum Rehab has equipped the unit with powerful new GLIDE motors for added power and torque.



ALSO AVAILABLE: *
Quantum 600 with XL
Package (400 lbs.
weight capacity)



QUANTUM® 610

The Quantum® 610 combines a compact frame with our patent pending Mid-Wheel 6® design platform to provide the ultimate choice for users who demand maximum stability, tight maneuverability and true all-around rehab performance. In addition, Quantum Rehab has equipped the unit with powerful new GLIDE motors for added power and torque.



* XL versions are available with Contoured Seat, Synergy® Seat, and TRU-Balance® Tilt only (dimensions up to 24" x 24").