



Spitfire II - 888WAL

SHOPRIDER®

The New Generation



Spitfire II - 888WAL

Shoprider's new rear wheel drive Spitfire II powerchair offers the best in style, comfort and performance. Enjoy a smooth ride while sitting in an elegant black leatherette seat with height adjustable armrests. The Spitfire II offers the ultimate in functionality and performance.

Built with the highest standards in quality and reliability

Since 1989 AMS HOMECARE INC. has built an enviable reputation by delivering high quality, reliable and comfortable self-care equipment addressing the needs of thousands of customers seeking to maintain their freedom and independence. We are dedicated to providing innovative, reasonably priced solutions specially developed with discriminating consumers in mind.

The Shoprider® brand represented exclusively in Canada by AMS HOMECARE INC. is recognized worldwide for superior design, technical sophistication as well as functionality in mobility products. Through warehouse-service centers conveniently located across the country, AMS HOMECARE INC. is the leader in the distribution of Shoprider® scooters and accessories.

In addition, we will be offering an expanding product line of AMS-branded products and accessories such as our current scooter canopy & cane holders and bed rails.

AMS HOMECARE INC. is a growing company committed to bringing to market the best products, devices and solutions to satisfy mobility and homecare needs.




Body style may differ from picture.


Features:

- Joystick Control
- Easy access on board charger
- Leatherette deluxe seat

Colours:



Blue



Red



Black

Specifications Spitfire II (model 888WAL):

Overall dimension: L* x W x H	mm (in)	1010 (40) x 680 (26.7) x 1168 (46)
Number, size of tires	Front	Front: 2 x 6"
	Rear	Rear: 2 x 14"
Weight capacity	kg (lb)	136.1 (300)
Total weight (with batteries, basket, seat)	kg (lb)	106 (233)
Battery capacity	_V_Ah x Pcs	12V50Ah x 2
Charger (on/off board)	V/A	4 Amp, on board
Output power of motor	hp	1000W/1.34 hp
Maximum speed	km/hr (mph)	9.9 (6.1)
Safe climbing angle	degree (%)	6 (10)
Max. driving distance with each full charge (depending on weight and terrain)		
	km (miles)	28 (17.5)
Turning radius	mm (inch)	660 (26)
Footrest		One piece with height, length and angle adjustment
Controller		Joystick, P & G programmable

All specification subject to change without prior notice.

Contact your AMS Homecare Dealer today:



AMS Homecare Inc. 1360 Cliveden Avenue, Delta, BC V3M 6K2 Canada
 Telephone: 604.273.5173 • Fax: 604.273.9312 • www.amshomecare.com

SHOPRIDER