



The Quickie Freestyle features an expansive array of options to ensure an optimal solution for any consumer!



Accessibility, freedom and style. The Quickie Freestyle™ has been thoughtfully designed to ensure a smooth ride. Our state-of-the-art power base provides unsurpassed obstacle climbing, intuitive maneuvering and a smooth ride.





With the Quickie Freestyle™ enjoy an active life with a sophisticated style. The Quickie Freestyle offers a complete range of options designed to meet your needs.



Quick-Change Shroud Colors
(4 colors provided w/ each chair)

Options:

- Readi-Align™ Casters
- Extended Range Battery Pack
- Walker Holder
- Cane/Crutch Holder
- O2 Holder
- Front Basket

Full Size F11 model



Features:

- Quick and responsive– High max speed (M11=5.5 mph) (F11=6.0 mph)
- Quick-change shroud available in red, blue, champagne & black (all colors included on all bases as standard)
- Compact and easy-to-maneuver design
- Overland™ All-Wheel Suspension & Readi-Align™ casters ensure obstacles are tackled with ease

FREESTYLE SPECIFICATIONS: M11 F11

Frame	M11	F11
Maximum Speed	5.5 mph	6.0 mph
Drive Wheels	10x3	14x3
Front Caster Size	6x2	6x2
Rear Suspension	S	S
Front Suspension	S	S
Overall Length w/footplate	38.5"	39.5"
Turning Radius	20"	22"
Obstacle Climbing	2.5"	3"
On-board Charger	4A	8A
Overall Width	22.5"	24.5"
HCPCS code	K0011	K0011

Control System	M11	F11
Integral	S	S
Qtronix Remote	0	0
Specialty Controls	0	0
Programmable	S	S

Battery Type	M11	F11
U1 (Gel)	0	
40 Ahr (AGM)	0	
22NF (Gel)		0
Gp 24 (Gel)		0

Max Weight Capacity	M11	F11
300 lbs	S	S
400 lbs		0

Seat Types	M11	F11
Rehab Seat Frame	0	0
High-back comfort	S	S

Comfort Seat Dimensions	M11	F11
16"x16"	S	S
18"x18"	S	S
20"x19"	S	S
22"x20"	S	S

Seat Frame Dimensions	M11	F11
Width	14"-24"	14"-24"
Depth	14"-20"	14"-20"
Back Height	15.5" 17.5" 19.5"	15.5" 17.5" 19.5"
Seat to Floor Height	15.5"-19.5"	17.5"-21.5"

Power Seating	M11	F11
Power C.G. Tilt (50°)	0	0
Smart Seat	na	na
Seat to Floor Height	17.5"-21.5"	19"-23"

Warranty	M11	F11
Frame	Lifetime	Lifetime
Integral controller	2 yrs	2 yrs
Remote Joystick	1 yr	1 yr
Remote Control Module	2 yrs	2 yrs
Motor	1 yr	1 yr
Gearbox	1 yr	1 yr

S– Standard 0– Optional upcharge