



## *Invacare*® | ***Toilet-Seat Risers***

Invacare Toilet Safety Products include a full range of products to make the bathroom a safer environment. The unique design, improved functions and outstanding quality make these models perfect compliments for home and care facility use. The safety you need and the quality you deserve, from the brand you trust.

Trusted Solutions,  
Passionate People



## Aquatec® Toilet Seat Risers



### A90

- High quality toilet seat riser.
- Raises the toilet seat by 4.1".
- Attaches safely to the toilet with two side clamps and four non-slip pads on base.
- High weight limit: 500 lbs.
- 2 year warranty.
- Item No. 1534999



### H430A

- For non-traditional or elongated bowls, the H430A is a secure, semi-permanent Raised Toilet Seat which offers superior safety and pleasing aesthetics for which the Aquatec brand is best known for.
- Item No.H430A



### H304 Finesse

- Toilet seat riser made of soft, waterproof latex free urethane
- Raises the toilet seat by 3.5".
- Fits on the toilet easily and safely, no assembly required.
- Large hygiene recess at the front for easy personal hygiene.
- Item No. 1575059



### Aquatec 90000

- Sturdy toilet seat riser.
- Available in two different seat heights: 2" & 4".
- Swing-away, ergonomic armrests for additional support when sitting down and standing up.
- Innovative clip system: The lid and toilet seat can be removed from the base; no tools are required.
- 2 year warranty.
- Item No. 1535002 - 2", 1535001 - 4"




### Aquatec 900

- High quality, all in one, multi height toilet seat riser.
- Individual adjustment. No tools required.
- Three seat heights - 2", 4" and 6".
- Three angles, lower at the front for ease of transition.
- Fold-way, ergonomic armrests for additional support when sitting down and standing up.
- Seat lid can be easily raised with the arm.
- 2 year warranty.
- Item No. 1535003


# Technical Data

Dimensions	A90	A900	A90000 - 2"	A90000 - 4"	Finesse	H430A
Seat Width	14"	14.1" (A)	14.4" (A)	14.4" (A)	15" (A)	14" (A)
Seat Depth	15.5"	16.7" (B)	15.4" (B)	15.4" (B)	16.7" (B)	15.2" (B)
Seat Height 1	4.1" (front)	2" (C)	2" (C)	4" (C)	3.5" (C)	2.3" (C)
Seat Height 2	4.5" (rear)	4" (D)				4" (D)
Seat Height 3		6" (E)				5.5" (E)
Total Width Armrests		23.2" (F)	21.2" (D)	21.2" (D)		22.2" (F)
Total Width No Armrests	15"	16.9" (G)	14.9" (E)	14.9" (E)	15.6"	15.5" (G)
Total Depth	16.1"	19.3" (H)	18.7" (F)	18.7" (F)	17" (F)	18.7" (H)
Total Height					5.3" (G)	13.8" (I)
Height Set to Armrest		7.9" (K)	7.9" (I)	7.9" (I)		7.8" (J)
Width Between Armrests		20" (L)	17.7" (J)	17.7" (J)		18.3" (K)
Width of Splash Guard		8.1" (M)				9" (L)
Depth of Splash Guard		10.2" (N)				9.1" (M)
Width of Opening	8.3"	8.9" (O)	8.1" (K)	8.1" (K)	8.8" (H)	10.1" (N)
Depth of Opening	10.6"	11.2" (P)	10.9" (L)	10.9" (L)	10.7" (I)	11.7" (O)
Weight Capacity	500 lbs.	265 lbs.	330 lbs.	330 lbs.	275 lbs.	265 lbs.


*Aquatec90*




*Finesse*




*Aquatec900*



*H430A*



*Aquatec90000*





### Molded Raised Toilet Seat - 1300RTS

- Attaches without tools and fits on most round bowl configurations.
- Integrated rear seat wings ensure secure attachment
- Contoured seat for internal rotation of hips.
- Simple to clean one piece molded design.

### Molded Clamp on Raised Toilet Seat - 1301RTS

- Clamp-on mechanism tightens easily with one hand.
- Comfortable contoured seat is made of easy to clean molded composite.
- Secures easily without tools to most round bowl configurations.

### Molded Clamp on Raised Toilet Seat with Arms - 1302RTS

- Adjustable padded locking arms provide additional stability
- Clamp-on mechanism tightens easily with one hand.
- Comfortable contoured seat is made of easy to clean molded composite.
- Secures easily without tools to most round bowl configurations.

## Technical Data

Model #	Dimensions with Arms	Dimensions Between Arms	Seat Dimensions	Seat-to-Arm Height	Arm-to-Floor Height	Product Weight	Weight Capacity
<b>1300RTS</b>	N/A	N/A	5" x 14.75" x 15.5" D	N/A	N/A	3 lbs. 6 oz.	300 lbs.
<b>1301RTS</b>	N/A	N/A	5" x 14.75" x 15.5" D	N/A	N/A	4.5 lbs.	300 lbs.
<b>1302RTS</b>	19.5"	N/A	5" x 14.75" x 15.5" D	10"	N/A	5.5 lbs.	300 lbs.



TRUSTED SOLUTIONS, PASSIONATE PEOPLE

Invacare Canada LP www.invacare.ca (800) 668-5324  
 570 Matheson Blvd. E., Unit 8, Mississauga, Ontario, L4Z 4G4 Canada. All rights reserved. Trademarks are identified by the symbols ™, ® and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. © 2017 Form No. 18-380C