INVACARE® TOP END®



CROSSFIRE™ ALL TERRAIN



WHEELCHAIR



The Invacare® Top End® Crossfire™ All Terrain Wheelchair is the monster truck of wheelchairs and is designed to leave some serious tracks. Do you like to hike, fish, camp, hunt, kayak, garden or go to the beach? You can do it all with the Crossfire All Terrain Wheelchair. Built of 6061 T6 Aluminum, it is equipped with enormous 25" x 2.75" knobby wheels and huge 8" x 2" front casters for extreme traction. The Crossfire All Terrain Wheelchair is engineered with mega ground clearance, plus forward and side-to-side stability so that you can conquer the toughest terrain.

Features

6061 T6 Aluminum frame
Adjustable center-of-gravity
Camber 0-9°
Lock down folding back
Adjustable back height
Adjustable tension seat upholstery
Adjustable tension back upholstery

24" or 25" high performance Sun® wheels
Hard anodized handrims
24" or 25" x 2.75" knobby tires
8" x 2" front casters
Adjustable height and angle footplate
Powder coat paint
250 lb. weight capacity

Model #: CRFAT



- 1. 24" or 25" x 2.75" knobby wheels and tires
- 2. 8" x 2" front casters
- 3. Engineered for front and side-to-side stability
- 4. Center-of-gravity adjustment
- 5. Anti-tip option
- 6. Optional push handles, wheel locks and anti-tip bars

Options

Spinergy® Spox wheels

Quad release axles

Titanium axles

Coated pushrims

Titanium handrims

Natural Fit® handrims

Fixed back

Push handles

Carbon fiber sideguards

Removeable sideguards

Folding sideguards

Wheel locks

Positioning straps

Alignment gauge

Footplate risers

2" cushion

Anti-tip bars

Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

Top End

www.topendwheelchair.com Sports and Recreation Products 4501 63rd Circle North Pinellas Park, Florida 33781-5914 (727) 522-8677 (800) 532-8677 Fax: (727) 522-1007

©2014

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols [™] and [®]. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Natural Fit is a registered trademark of Three Rivers Out Front. Spinergy is a registered trademark of Spinergy, Inc. Specifications are subject to change without notification. Form No. 12-504 Rev. 02/14 140182

