

 **TOP END**[®]



Invacare[®] Top End[®] XLT Handcycle

If recreation is your goal, the **Invacare[®] Top End[®] XLT Handcycle** is great fun and features low maintenance 7 speed hubs with reverse hands-on braking. Just get out there and go! All Invacare Top End recreation handcycles are equipped with a safety flag, full composite chain guard, cushion and seat/chest restraint, as standard, for extra value, comfort and safety.

The **Invacare® Top End® XLT Handcycle** is designed with the recreational handcycling enthusiast in mind. Equipped with seven speeds, the XLT handcycle will easily handle most terrain, whether it's riding around town or taking a spin in the country. The lightweight, oval tube aluminum XLT frame with chrome-moly steel fork features simple geometry for maximum stiffness and easy cornering. The seat and back are totally adjustable for a comfortable fit. The narrow rear-end enhances aerodynamics and the 15° of camber increases stability. If you want a great way to exercise, cross-train or just have fun, choose the XLT handcycle.

Model #: XLT



1. Twist Shifter
2. Quad Twist shifting option
3. QuadGrip handpedals with 3" width extension option



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End
www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(800) 532-8677
Fax: (727) 522-1007

©2018
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. Shimano
is a registered trademark of
Shimano Inc. Specifications
are subject to change
without notification.
Form No. 18-323
181125

Features

- Choice of 15" wide or 17" wide seat widths
- 6061 Aluminum frame/chrome-moly fork
- Adjustable seat, back angle and crank height
- Full chain guard, chain tension idler and safety flag
- Adjustable footrest system with safety straps
- Seat and back cushions with chest restraint
- Shimano® 7 speed hub with reversing drum brake
- 175 mm crankset with one chain ring
- Quad Twist shifter option helps those with hand impairments
- Frictionless ovalized ergonomic hand pedals
- Parking brake 26" High Performance wheels
- High pressure clincher tires
- Quick release axles
- 15° camber
- Road crown compensator
- Reflectors
- 250 lb. weight capacity
- Many options and aftermarket accessories are available
- Mountain Drive
- Quad Grips
- Tri-Pin Pedals

TOP END®