

INVACARE® MATRX®
SEATING SERIES
 PRODUCT GUIDE



Invacare® Matrx® Seating and Positioning is an established global leader, committed to product innovation, high quality standards and superior product performance. As the world of seating evolves, the Invacare Matrx Seating Series continues to exceed customer expectations.

Call to order

1.800.668.5324
 invacare.ca

1.866.748.7943
 motionconcepts.com



matrx®
 SEATING SERIES

TABLE OF CONTENTS

Pediatric

Cushions

PS	4
PSP	4
Kid*ab*ra Vi	4

Backs

Mini Back	5
Mini Fixed Laterals	5
Mini Swing-away Laterals	5
Mini Chest Harness Interface	5
MiniEXT	5

Accessories

Mini Elan Headrest	4
--------------------	---

Standard

Cushions

Libra™ / Libra™ Adjust	6
Libra™ Fit Kit	6
PS	7
PSP	7
PSVF	7
Vi	7
Stabilite™	7
Flovair™ Gentle/Max Contour	7
Options	8
General Use	8

Backs

Elite™	9
Elite™ Deep	9
Elite™ TR	9
PB	9
PB Deep	9
Hardware	10
Options	10
Accessories	11
General Use	8

High Performance

MX2 Back	12
----------	----

Bariatric

Cushions

PS Heavy-Duty	14
Vi Heavy-Duty	14
Libra™ Heavy-Duty/ Libra™ Adjust Heavy-Duty	14
Libra™ Fit Kit	6
Flovair™ Heavy-Duty Gentle/Max Contour	14
Stabilite™ Heavy-Duty	14

Backs

Elite™ Heavy-Duty	15
Elite™ Deep Heavy-Duty	15
Elite™ TR Heavy-Duty	15
PB Heavy-Duty	15
PB Deep Heavy-Duty	15

Accessories

Elan Headrest	16
Arm Pads/Straps/Lap Belts/Trays	17

Quick Product Reference

Quick Product Reference	18/19
-------------------------	-------

INVACARE® MATRX® MINI CUSHIONS



PS Pediatric Cushion

ADP #: SE0000131 (10" wide & smaller)
SES1M0065 (11"-12" wide)
SES1M0060 (13" wide & up)

PS Pediatric Cushion is an exceptionally lightweight cushion designed to provide superior positioning, stability and comfort for our smallest clients. The design features an anatomically contoured shape with a unique waffled ischial relief recess and an innovative reversible outer cover.

Width: 8" - 14"

Depth: *10", 12", 14", 16"

* available for 8"-10" W cushions only

Weight Capacity: 150 lb.

Fabric Options: Kid*ab*ra Meshtex, Blizzard Blue, Pink Lollipop and Black Startex (shown)

Options: Custom Modifications, Cushion Rigidizer, Pelvic Obliquity Kit

■ Available only from Motion Concepts



PSP Pediatric Cushion

ADP #: Pending

The Matrix PSP (Posture Seat Polymer) cushion combines the stability and support of the PS foam base with a polymer overlay for extra protection and shear reduction. An ultra-lightweight cushion providing zero maintenance, worry free skin protection and positioning.

Width: 14"

Depth: 12", 14", 16"

Weight Capacity: 150 lb.

Fabric Options: Kid*ab*ra Meshtex, Blizzard Blue, Pink Lollipop and Black Startex (shown)

Options: Custom Modifications, Cushion Rigidizer, Pelvic Obliquity Kit

■ Available only from Motion Concepts



Kid*ab*ra Vi Cushion

ADP #: SES1M0015

Kid*ab*ra Vi Cushion is designed to provide skin protection and positioning for individuals at moderate to high risk of skin breakdown. The unique three-part foam construction targets skin protection where it is needed most, and features excellent stability, support and comfort.

Width: 12", 14"

Depth: 12", 14"

Weight Capacity: 150 lb.

Cover Options: Kid*ab*ra Meshtex (shown), Black Meshtex

Options: Custom Modifications, Drop Base, Cushion Rigidizer, Pelvic Obliquity Kit

INVACARE® MATRX® ELAN HEADREST



Mini Elan Headrest

ADP #: SEHND0015 (Hardware) +
SE0001320 (Standard Pad w/Startex Cover)
SE0001321 (Standard Pad w/Infection Control Cover)
SE0001324 (4 Point Pad w/Startex Cover)
SE0001325 (4 Point Pad w/Infection Control Cover)

Compact, lightweight and highly adjustable

Adjustable Headrest Hardware:

Six points of adjustment through pivot ball and multi-disk locking plates; includes tool-less adjustment levers. Vertical Post 11" High.

Pads: Pediatric pad options available, Standard Startex or Infection Control covers; aluminum plates can be formed to shape. Other pad options available on Page 16



6" x 3.5" Standard Pad



9" x 7" 4 Point Pad

INVACARE® MATRX® MINI BACK



Mini Back

ADP #: SE0001213

Mini Back features Mini Set hardware. The Mini Back is contoured to fit children of all shapes and sizes, combining exceptional support in an appealing low-profile package.

Width: 10", 11", 12", 13", 14"

Height: 7", 10", 12", 14"

Weight Capacity: 150 lb.

Contour Depth: 3"

Fabric Options: Kid*ab*ra Meshtex (shown); Blizzard Blue, Pink Lollipop and Black Startex

Standard Hardware: Mini-Set Mounting Hardware (Quick-release latch for easy removal/portability; Fits 1", 7/8" and 3/4" diameter canes; allows back angle/ depth adjustment with child seated in wheelchair)

Options: Swing-away Laterals, Fixed Laterals, Chest Harness Interface, MiniEXTHardware Interface

■ Available only from Motion Concepts



Mini Fixed Laterals

ADP #: SELD0005

Designed especially for the Mini Back. Installs inside back shell; Does not compromise depth/angle adjustment.

Pad Size: 4.5" L x 3" H

Fabric Options: Kid*ab*ra Meshtex (shown); Blizzard Blue, Pink Lollipop and Black Startex

■ Available only from Motion Concepts



Mini Swing-Away Laterals

ADP #: SEICF202L (Hardware) + SELCF102L (Pad)

Designed for the Invacare Matrix Mini Back. Installs inside back shell; Does not compromise depth/angle adjustment. Lock/release latch allows pad to swing away for side transfers.

Pad Size: 4.25" L x 3.25" H, 5.25" L x 3.5" H

Fabric Options: Black Startex

■ Available only from Motion Concepts

Pediatric Cushions/Backs Fabric Options



Pink Lollipop Startex



Blizzard Blue Startex



Kid*ab*ra Meshtex



Black Startex



Chest Harness Interface (HIMINI)

Height adjustable hardware allows attachment of OEM-compliant abdominal supports to the full range of Mini backs.

■ Available only from Motion Concepts



MiniEXT Hardware

Provides pin mounting location +/- 3.5" for clearance of obstacles on back canes.

■ Available only from Motion Concepts



Mini Set Hardware

Standard on all Pediatric Elite and Posture Backs.

Provides 40° Angle, 4" depth and even rotation adjustment.

2015 AWARD
HARDING
 INNOVATIVE
 PRODUCT
CMC15
 SEATING, MOBILITY AND MORE

INVACARE® MATRX® CUSHIONS

Libra™ Cushions

ADP #: SE0001212

Libra™ Cushion is designed with the highest level of skin protection, positioning and adjustability for unique user needs. It is specially contoured to provide balance and stability to the pelvis, optimizing body function. Recommended for all users, including those at highest risk of skin breakdown.

Width: 14", 15", 16", 17", 18", 19", 20"

Depth: *12", *14", 16", 18", 20", 22"

* available for 14" W cushions only

Weight Capacity: 300 lb.

Options: Custom Modifications, Drop Base, Cushion Rigidizer, Pelvic Obliquity Kit, Fit Kit Accessories
 Custom options available through Motion Concepts only


Cover Options: Startex™ (shown), Spacer Mesh, Infection Control (wipeable)



Libra Fit Kit Accessories

ADP #: SESND1040/SESND1045

Designed to enhance immersion, correct or accommodate orthopedic issues



new



Libra Adjust Cushions

ADP #: SE0001212 + SESND1045

Libra™ Adjust Cushion includes 2 Fluid Supplement Packs to create deep contour pelvic well. Recommended for taller individuals or those with a high degree of muscle wasting.

Width: 14", 15", 16", 17", 18", 19", 20"

Depth: *12", *14", 16", 18", 20", 22"


* available for 14" W cushions only

Weight Capacity: 300 lb.

Options: Custom Modifications, Drop Base, Cushion Rigidizer, Pelvic Obliquity Kit, Fit Kit Accessories
 Custom options available through Motion Concepts only

Cover Options: Startex™ (shown), Spacer Mesh, Infection Control (wipeable)

■ Available only from Motion Concepts





PS Cushion

ADP #: SES1M0060

PS Cushion is an exceptionally lightweight cushion designed to provide superior positioning, stability and comfort. The design features an anatomically contoured shape with a unique waffled ischial relief recess and an innovative reversible outer cover. It is a great choice for a very wide range of individual needs.

Width: *14", 16", 17", 18", 20"

* available for 16" D cushions only

Depth: 16", 18", 20", 22"

Weight Capacity: 300 lb.

Options: Custom Modifications, Cushion Rigidizer, Pelvic Obliquity Kit



PSP Cushion

ADP #: Pending

PSP (Posture Seat Polymer) cushion combines the stability and support of the PS foam base with a polymer overlay for extra protection and shear reduction. An ultra-lightweight cushion providing zero maintenance, worry free skin protection and positioning.

Width: *14", 16", 17", 18", 20"

* available for 16" D cushions only

Depth: 16", 18", 20", 22"

Weight Capacity: 300 lb.

Options: Custom Modifications, Cushion Rigidizer, Pelvic Obliquity Kit

■ Available only from Motion Concepts



PSVF Cushion

ADP #: SES1M0070

The unique three part foam construction combination of PSVF cushion provides maximum immersion and comfort without sacrificing long lasting durability and support. Exceptionally comfortable, a popular choice for geriatric or sensitive users.

Width: 16", 18", 20"

Depth: 16", 18", 20"

Weight Capacity: 300 lb.

Options: Custom Modifications, Cushion Rigidizer, Pelvic Obliquity Kit



Vi Cushion

ADP #: SES1M0005/ SES1M0010

Vi Cushion is designed to provide skin protection and positioning for individuals at moderate to high risk of skin breakdown. The unique three-part foam construction targets skin protection where it is needed most, and features excellent stability, support and comfort.

Width: 14", 16", 17", 18", 20"

Depth: *15", *16", 17", 18", 20", 22"

* available for 14" W cushions only

Weight Capacity: 300 lb.

Options: Custom Modifications, Drop Base, Cushion Rigidizer, Pelvic Obliquity Kit



Invacare® Stabilite™ Cushion

ADP #: SE0000966

Stabilite Cushion comes standard with an Invacare ThinAir Bladder to assist in redistributing peak pressure away from bony prominences and help manage shear forces. Contoured, slim cushion insert helps eliminate sling upholstery "hammock effect."

Width: 14" - 20"

Depth: 14" - 20"

Weight Capacity: 250 lb.

Options: Custom Modifications, Cushion Rigidizer, Removable Leg Wedges, Pelvic Obliquity Kit

■ Available only from Invacare



Invacare® Flovair™ Gentle/Max Contour Cushion

ADP #: SE0000967

Flovair Gentle Contour cushion captures the stability of foam, the pressure redistribution of air and the shear management of fluid all in one hybrid cushion. Zippered inner cover easily wipes clean and allows access for cushion modifications.

Width: 15" - 20"

Depth: 15" - 20"

Weight Capacity: 300 lb.

Options: Custom Modifications, Cushion Rigidizer, Removable Leg Wedges, Pelvic Obliquity Kit

■ Available only from Invacare

INVACARE® MATRIX® OPTIONS



Cushion & Cover Modifications

Custom modifications are available to accommodate a variety of orthopedic conditions. Standard modifications include seat width and depth reductions, as well as leg length discrepancies. (Additional lead time required)



Pelvic Obliquity Kit (2 pcs)

Designed to accommodate up to 1° obliquities. Can be readily modified and installed onto the following cushions:

Vi and Kid*ab*ra Vi Cushions; PS Cushions; PSP Cushions; Flovair Cushions; Stabilite Cushions



Cushion Rigidizers

Lightweight plastic seat insert fits between seat rails for level support and eliminates hammock effect of sling seat upholstery. Securely fitted to bottom of foam cushion inside the cushion cover.

INVACARE® MATRX® BACKS

Standard



Elite™ Back

ADP #: SEB1M0025

Elite Back utilizes a contoured aluminum back shell with scapular cutout. MeshTex cover and foam overlap back shell for comfort and protection.

Width: 14" - 20"

Height: 10", 12", 14", 16", 18", 20"

Weight Capacity: 300 lb.

Contour Depth: 3"

Standard Hardware: Easy Set

Hardware Options: Extended Range Hardware, Split Hardware, Fixed Space Saver*, Heavy-Duty Easy Set



Elite™ Deep Back

ADP #: SEB1M0026

Elite Deep Back designed with extra contour depth for those requiring more midline support.

Width: 14" - 20"

Height: 10", 12", 14", 16", 18", 20"

Weight Capacity: 300 lb.

Contour Depth: 6"

Standard Hardware: Easy Set

Hardware Options: Extended Range Hardware, Split Hardware, Fixed Space Saver*, Heavy-Duty Easy Set



Elite™ TR Back

ADP #: SE0001214

Elite TR Back is ideal for Tilt or Tilt/Recline applications. Includes SuperSoft foam for extra immersion & comfort.

Width: 15", 16", 17", 18", 19", 20"

Height: 16", 18", 20"

Weight Capacity: 300 lb.

Contour Depth: 3"

Made with SuperSoft Foam: (Dual layer foam pad with 1" super-soft foam laminated to 1" firm HR foam)

Standard Hardware: Easy Set

Hardware Options: Extended Range Hardware, Split Hardware, Fixed Space Saver*, Heavy-Duty Easy Set



PB Back

ADP #: SEB1M0005

Elite PB Back utilizes an ABS Shell with molded HR foam. Gently contoured shape with full wraparound upholstery.

Width: 14", 16", 17", 18", 20"

Height: 12", 16", 20"

Weight Capacity: 300 lb.

Contour Depth: 5"

Standard Hardware: Easy Set

Hardware Options: Extended Range Hardware, Split Hardware, Fixed Space Saver*, Heavy-Duty Easy Set



PB Deep Back

ADP #: SEB1M0015

PB Deep Back offers extra deep contour for midline support and stability.

Width: 14", 16", 17", 18", 20"

Height: 12", 16", 20"

Weight Capacity: 300 lb.

Contour Depth: 7"

Standard Hardware: Easy Set

Hardware Options: Extended Range Hardware, Split Hardware, Fixed Space Saver*, Heavy-Duty Easy Set



Easy Set Hardware

Standard on all Elite and Posture Backs.

Provides 40° Angle,

4" depth and

even rotation adjustment.

For other hardware

options refer to Page 10

**Note: This back can be installed on the Invacare® Formula™ CG Tilt/Recline, and Motion Concepts TRx Power Positioning Systems using FSSH*

INVACARE® MATRX® HARDWARE



Mounting Hardware Options

Easy Set (PBMH)

300 lb. weight capacity. Quick release design for easy removal/portability; Allows back angle/depth adjustment with client seated in wheelchair; Compatible with all models of Elite & PB backs.

Split (SMH)

500 lb. weight capacity. Allows clearance of back cane rigidizer bars or other obstacles.

Heavy-Duty Easy Set (HDMH)

500 lb. weight capacity. Included with all Heavy-Duty Back models; Recommended for users with high tone or heavy intended use.

Fixed Space Saver (FSSH)

500 lb. weight capacity. Recommended for power tilt/recline or tighter width applications. Permanent mount requires tool for removal.

Mounting Pins

Standard/Long Mounting Pins: Included with PB and Elite Backs; allows PB Back to fit 1" narrower to 1" wider and Elite backs the same width or 1" wider wheelchairs.

Optional Extended Mounting Pins (EMP): Allow interface to wheelchairs up to 2" wider. 300 lb. weight capacity maximum.

INVACARE® MATRX® OPTIONS



Super-soft Foam Upgrade

Provides a 1" layer of super-soft foam laminated onto a 1" layer of firmer high resilience foam for increased immersion and added comfort. Available for the following models/sizes:

Invacare Matrix Elite & Elite Deep Backs

Width: 16", 18", 20"

Height: 16", 18", 20"



Visco Foam Upgrade

Provides a ¾" layer of visco foam laminated onto a 1" layer of high resilience foam. for increased immersion and added comfort. Visco upgrade is available on all Invacare Matrix PB back models.

(Note: option not available with molded foam)



Waterproof Fabric

Wipeable Slip Cover

Additional cover adds a layer of protection for the back cover and foam. Easily removed for cleaning and ideal for incontinent users and institutions. Available for all Elite and PB Back models.

INVACARE® MATRIX® BACK ACCESSORIES



Fixed Laterals

ADP #: SELD0005

Designed for Invacare Matrix Elite & Elite Heavy-Duty Backs. Mounting: Installs inside back shell; Does not compromise depth/angle adjustment.

Options: Standard/ Pediatric/Offset (Left or Right)

Sizes:

- Standard 7.5" L x 5" H
- Pediatric 5.5" L x 3" H
- Offset 7.5" L x 5" H



Swing Away Laterals

ADP #: SEICF202L (Hardware) + SELCF102L (Pad)

Designed for Invacare Matrix Elite and Elite Heavy-Duty Backs. Mounting: Installs inside back shell; Does not compromise depth/angle adjustment.

Lock/release latch allows pad to swing away for side transfers.

Options: Standard or Offset (Left or Right)

Pad Sizes:

- 4.25" L x 3.25" H
- 5.25" L x 3.5" H
- 5.75" L x 3.75" H
- 6.5" L x 4.25" H



Chest Harness Interface for Elite Backs (CHIEB)

Width, height and angle adjustable chest harness interface for precise orientation of shoulder straps.



Chest Harness Interface for PB Backs (CHIPB)

Width adjustable chest harness interface for precise orientation of shoulder straps.



Privacy Flaps (PFLAP)

Available for all models of Invacare Matrix PB, Elite (shown), and the MX1 backs.

Standard Sizes: XS up to XXL, (fit back widths 13" up to 24")
Note: Flaps over 24" wide are available as special order only

Standard



SEMPER
PARATUS

SCHWALBE

K-SERIES
attract

2017 AWARD
HARDING
 INNOVATIVE
 PRODUCT
CMC17
 SEATING, MOBILITY AND MORE



INVACARE® MATRX® MX2 BACK



Matrx MX2 Back

- ADP #: SE0001460** (all 6" High)
- SE0001461** (all 9" High)
- SE0001462** (all 12" High)
- SE0001463** (all 16" High)

MX2 Back is designed for active users who demand the best. Carbon fiber specially engineered to provide unmatched stability for peak performance. Contoured to provide hip/trunk stability and enhance upper extremity mobility and function. Very light weight and aesthetically appealing.

Width: 14", 15", 16", 17", 18"

Height: 6", 9", 12", 16"

Weight Capacity: 250 lb.

Contour Depth: 3"

Hardware Options: Fixed Mount & Quick Release

■ Available only from Motion Concepts

Solid Carbon Fiber Back Shell construction designed to provide optimal flex and suspension. Shaped for perfect fit and maximum stability. Designed with cutouts for breathability and slots for easy height adjustment.



Fixed Mount



Quick Release

Hardware Options

Angle, height, depth, width, rotation adjustable for precise orientation. Width adjustable to fit chairs 1" narrower or wider than back model widths.

Fixed Mount – Lightweight, carbon fiber construction, maximum performance for rigid chairs

Quick Release – Easy removal for portability



MX2 Deep Laterals

Provide extra mid-line support and stability.

Pad Size: 4.5" L x 3" H

■ Available only from Motion Concepts



Privacy Flaps (PFLAP)

Available for all MX2 back sizes

High Performance

INVACARE® MATRX® HEAVY-DUTY CUSHIONS



PS Heavy-Duty

ADP #: SES1M0066

PS Cushion is designed to provide superior positioning, stability and comfort. Features an anatomically contoured shape with a unique waffled ischial relief recess and an innovative reversible outer cover.

Width: 21" - 30"

Depth: 18" - 30"

Weight Capacity: 600 lb.

Options: Custom Modifications, Pelvic Obliquity Kit, Cushion Rigidizer



Vi Heavy-Duty Cushion

ADP #: SES1M0011

Vi Cushion is designed to provide skin protection and positioning for individuals at moderate to high risk of skin breakdown. The unique three-part foam construction targets skin protection where it is needed most, and features excellent stability, support and comfort.

Width: 21" - 30"

Depth: 21" - 26"

Weight Capacity: 600 lb.

Options: Custom Modifications, Pelvic Obliquity Kit



Libra Heavy-Duty Cushion

ADP #: SE0001237

Libra™ Cushion is specially contoured to provide balance and stability to the pelvis, optimizing body function. Recommended for all users, including those at highest risk of skin breakdown

Width: 21", 22", 23", 24"

Depth: 18", 20", 22"

Weight Capacity: 500 lb.

Cover Options: Startex™ (shown), Spacer Mesh, Infection Control (wipeable)

Options: Fit Kit Accessories Custom options available through Motion Concepts only



Invacare® Flovair™ Heavy-Duty Gentle/Max Contour Cushion

ADP #: SESND1060

Flovair Gentle Contour cushion captures the stability of foam, the pressure redistribution of air and the shear management of fluid all in one hybrid cushion.

Width: 21" - 24"

Depth: 16" - 22"

Weight Capacity: 500 lb.

Options: Custom Modifications, Cushion Rigidizer, Removable Leg Wedges, Pelvic Obliquity Kit

Available only from Invacare



Invacare® Stabilite™ Heavy-Duty Cushion

ADP #: SESND1060

Stabilite Cushion comes standard with an Invacare ThinAir Bladder to assist in redistributing peak pressure away from bony prominences and help manage shear forces.

Width: 21" - 24"

Depth: 16" - 22"

Weight Capacity: 350 lb.

Options: Custom Modifications, Cushion Rigidizer, Removable Leg Wedges, Pelvic Obliquity Kit

Available only from Invacare



Libra™ Adjust Heavy-Duty Cushion

ADP #: SE0001237+SESND1045

Libra™ Adjust Cushion features 2 Fluid Supplement Packs to create deep contour pelvic well. Recommended for taller individuals or those with a high degree of muscle wasting.

Width: 21", 22", 23", 24"

Depth: 18", 20", 22"

Weight Capacity: 500 lb.

Cover Options: Startex™ (shown), Spacer Mesh, Infection Control (wipeable)

Options: Fit Kit Accessories Custom options

Available only from Motion Concepts

INVACARE® MATRIX® HEAVY-DUTY BACKS

Bariatric



Elite™ Heavy-Duty Back

ADP #: SE0001051

Elite HD Back utilizes reinforced aluminum back shell for Bariatric users.

Width: 21" - 30"

Height: 16", 20"

Weight Capacity: 500 lb.

Extra Wide PB Heavy-Duty Backs are available up to 30" W x 16" or 20" H

Contour Depth: 3"

Standard Hardware: Heavy-Duty Easy Set

Hardware Options: Fixed Space Saver*



Elite™ Deep Heavy-Duty Back

ADP #: SEB1M0040

Elite Deep HD Back is designed with extra contour depth for those requiring midline support.

Width: 21" - 24"

Height: 16", 18", 20"

Weight Capacity: 500 lb.

Contour Depth: 6"

Standard Hardware: Heavy-Duty Easy Set

Hardware Options: Fixed Space Saver*



Elite™ TR Heavy-Duty Back

ADP #: SE0001051+

SEBND0020 (Additional pressure relief)

Elite TR HD Back is ideal for Tilt or Tilt/Recline applications. Includes SuperSoft foam for extra immersion & comfort.

Width: 21" - 24"

Height: 20"

Weight Capacity: 500 lb.

Contour Depth: 3"

Made with SuperSoft Foam: (Dual density, 1" firm + 1" soft HR foam combination)

Standard Hardware: Heavy-Duty Easy Set

Hardware Options: Fixed Space Saver*



PB Heavy-Duty Back

ADP #: SEB1M0010

PB HD Back utilizes contoured aluminum shell for Bariatric users. Full wraparound upholstery.

Width: 20", 22", 24"

Height: 16", 20"

Weight Capacity: 500 lb.

Extra Wide PB Heavy-Duty Backs are available up to 30" W x 16" or 20" H

Contour Depth: 5"

Standard Hardware: Heavy-Duty Easy Set

Hardware Options: Fixed Space Saver*



PB Deep Heavy-Duty Back

ADP #: SEB1M0020

PB Deep HD Back offers extra deep contour for midline support and stability.

Width: 20", 22", 24"

Height: 16", 20"

Weight Capacity: 500 lb.

Extra Wide PB Deep Heavy-Duty Backs are available up to 30" W x 16" or 20" H

Contour Depth: 7"

Standard Hardware: Heavy-Duty Easy Set

Hardware Options: Fixed Space Saver*



Heavy-Duty Easy Set Hardware

Standard on all Heavy-Duty Elite and Posture Backs.

Provides 40° Angle,

4" depth and

even rotation adjustment.


For other hardware

options refer to Page 10

*Note: This back can be installed on the Invacare® Formula™ CG Tilt/Recline, and Motion Concepts TRx Power Positioning Systems using FSSH



INVACARE® MATRX® ELAN® HEADREST



Elan Headrest Hardware

ADP #: SEHND0015 (Hardware)

Compact, lightweight and highly adjustable

Elan Hardware: Six points of adjustment through pivot ball and multi-disk locking plates, available in standard -15" post or mini -11" post; includes tool-less adjustment levers

Pads: 6 Pad options available; aluminum plates can be formed to shape

Pad Covers:

Standard Startex or Infection Control covers available

■ Available only from Motion Concepts



Co-molded Foam



Tool-less Adjustment levers



Ergonomic Knob



Standard Pad

ADP #: SE0001320 (Startex Cover)
SE0001321 (Infection Control Cover)

Pad Sizes:

- 6" Pad: 6" x 3.5"
- 10" Pad: 10" x 5"
- 14" Pad: 14" x 5"



Occipital Pad

ADP #: SE0001322 (Startex Cover)
SE0001323 (Infection Control Cover)

Pad Sizes:

- 9" Pad: 9" x 7"
- 12" Pad: 12" x 8"



4 Point Pad

ADP #: SE0001324 (Startex Cover)
SE0001325 (Infection Control Cover)

Pad Size:

- Small Pad: 9" x 7"
- Large Pad: 11" x 10"

INVACARE® MATRIX® ACCESSORIES



Ergonomic Arm Trough

ADP #: SEAND0005 (Arm Troughs -max 2)
SEAND0015 (Hand Pad -max 2)

Molded, ergonomic shape fits the natural arm curve; provides users flexibility of various arm positions Available for left or right side and is 13" x 4.5"; hand pad is optional add-on and is also available for left or right side and is 5" x 4.5"

Multi-Axis Armrest Mount

ADP #: SEAND0025 (max 2)

Allows infinite adjustability of angle and rotation of arm pad or Multi-position arm pads or Arm Trough.



Modular Arm Pads

ADP #: SEAND0030 (Visco Foam -max 2)
SEAND0010 (Integral Foam - max 2)

Durable composite frame features 3 slots for infinite fore and aft adjustability as well as offset mounting.

Available for left or right side and in 2 sizes:

Full length: 14" x 4"

Desk length: 10" x 4"



Gel Arm Pads

ADP #: SEAND0030 (max 2)

Available for left or right side and in 3 sizes:

Narrow: 12" x 2"

Wide: 12" x 3"

Long: 14" x 3.5"



Multi-Position Arm Pad

ADP #: SEAND0010 (max 2)

Available for left or right side and in 2 sizes:

Full length: 14" x 3"

Desk length: 10" x 3"



Bodypoint Padded Lap Belts

ADP #: SERND0030 - 2 Point
SERND0020 - 4 Point

Available in 2 and 4 Point; Central pull, push release bottom; fits hip widths 13" - 19" W



Bodypoint Aeromesh Calf Straps/Panels

ADP #: SEFCF107L

Available in 4 sizes:

3.5" H and fits chairs 15" - 18" W

3.5" H and fits chairs 18" - 22" W

9" H and fits chairs 15" - 18" W

9" H and fits chairs 18" - 22" W



Invacare® Drop Base

Provides stable support base for the cushion. Three sizes fit wheelchairs from 14" - 22" W

■ Available only from Invacare

Cushions	Stabilite	Stabilite Heavy-Duty	PS Pediatric	PS	PS Heavy-Duty	PSP	PSP	PSVF
	STANDARD	BARIATRIC	PEDIATRIC	STANDARD	BARIATRIC	PEDIATRIC	STANDARD	STANDARD
Model	ITSGwd	ITSGwd	PSwwdd	PSwwdd	PSwwdd	PSwwdd	PSwwdd	PSwwdd
Width	14" - 20"	21" - 24"	08" - 14"	16", 18", 20"	21" - 30"	14"	14", 16", 17", 18", 20"	16", 18", 20"
Depth	14" - 22"	16" - 22"	10", 12", 14", 16"	16", 18", 20", 22"	18" - 30"	12", 14"	16", 18", 20", 22"	16", 18", 20"
Weight Capacity	250 lb.	350 lb.	150 lb.	300 lb.	600 lb.	150 lb.	300 lb.	300 lb.
Std. Cover	Meshtex	Meshtex	Startex (Reversible)	Startex (Reversible)	Startex (Reversible)	Startex (Reversible)	Startex (Reversible)	Startex (Reversible)
Cover Options			Kid*a*bra, Pink Lollipop, Blizzard Blue, Black Startex			Kid*a*bra, Pink Lollipop, Blizzard Blue, Black Startex		
ADP Codes	SE0000966	SESND1060	SE0000131 SES1M0065 SES1M0060	SES1M0060	SES1M0066	PENDING	PENDING	SES1M0070
Available from	INVACARE	INVACARE	INVACARE	INVACARE	INVACARE	MOTION CONCEPTS	MOTION CONCEPTS	INVACARE
			MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS			MOTION CONCEPTS

Backs	Mini	Elite	Elite Heavy-Duty	Elite TR/ Elite SuperSoft®	Elite TR Heavy-Duty	Elite Deep	Elite Deep Heavy-Duty	PB
	PEDIATRIC	STANDARD	BARIATRIC	STANDARD	BARIATRIC	STANDARD	BARIATRIC	STANDARD
Model	MINIwwhh	PBEwwhh	PBEwwhh-HD	PBEwwhh-TR	PBEwwhh-TRHD	EDBwwhh	EDBwwhh-HD	MPBwwhh
Width	10" - 14"	14" - 20"	21" - 24" ; 25" - 30" (Extra Wide)	15" - 20" 16", 18", 20" *	21" - 24"	14" - 20"	21" - 24"	14", 16", 17", 18", 20"
Height	7", 10", 12", 14"	10", 12", 14" 16", 18", 20"	16", 20"	16", 18", 20"	20"	10", 12", 14" 16", 18", 20"	16", 18", 20"	12", 16", 20"
Contour Depth	3"	3"	3"	3"	3"	6"	6"	5"
Weight Capacity	150 lb.	300 lb.	500 lb.	300 lb.	500 lb.	300 lb.	500 lb.	300 lb.
Std. Cover	Kid*a*bra	Meshtex	Meshtex	Meshtex	Meshtex	Meshtex	Meshtex	Meshtex
Cover Options	Pink Lollipop, Blizzard Blue	Slip Cover (wipeable)**		Slip Cover (wipeable)**		Slip Cover (wipeable)**		
Fits onto Back Canes:	Fits up to 1" wider	Fits up to 1" wider	Fits up to 1" wider	Fits up to 1" wider	Fits up to 1" wider	Fits up to 1" wider	Fits up to 1" wider	Fits 1" narrower up to 1" wider
ADP Codes	SE0001213	SEB1M0025	SE0001051	SE0001214	SE0001051	SEB1M0026	SEB1M00401	SEB1M0005
Available from	MOTION CONCEPTS	INVACARE	INVACARE	INVACARE	INVACARE	INVACARE	INVACARE	INVACARE
		MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS

*This information is not intended to be, nor should it be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

Flovair Gentle/Max Contour	Flovair Heavy-Duty Gentle/Max Contour	Kid*ab*ra-Vi	Vi	Vi Heavy-Duty	Libra Libra Adjust	Libra Heavy-Duty Libra Adjust Heavy-Duty
STANDARD	BARIATRIC	PEDIATRIC	STANDARD	BARIATRIC	STANDARD	BARIATRIC
ITFGwd /ITFMwd	ITFGwd /ITFMwd	KAwwdd-Vi	MAwwdd-Vi	MAwwdd-Vi	LCwwdd	LCwwdd
15" - 20"	21" - 24"	12", 14"	14", 16", 17", 18", 20"	21" - 30"	14" - 20"	16" - 24"
15" - 20"	16" - 22"	12", 14"	15", 16", 17", 18", 20", 22"	18" - 30"	12", 14", 16", 18", 20", 22"	18", 20", 22"
300 lb.	500 lb.	150 lb.	300 lb.	600 lb.	300 lb.	500 lb.
Meshtex	Meshtex	Kid*a*bra	Meshtex	Meshtex	Startex	Startex
		Meshtex			Spacetex (-SP), Infection Control (-IC)*	Spacetex (-SP), Infection Control (-IC)*
SE0000967	SESND1060	SES1M0015	SES1M0005 SES1M0010	SES1M0011	SE0001212 SE0001212+ SESND1045	SE0001237 SE0001237+ SESND1045
INVACARE	INVACARE	INVACARE	INVACARE	INVACARE	INVACARE	INVACARE
		MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS

PB Heavy-Duty	PB Deep	PB Deep Heavy-Duty	MX2 Fixed
BARIATRIC	STANDARD	BARIATRIC	STANDARD
HDBwwhh	PBDwwhh	DHDwwhh	IMX2wwhh-CFH
20" - 30" †	14", 16", 17", 18", 20"	20" - 30" †	14" - 18"
16", 20"	12", 16", 20"	16", 20"	6", 9", 12", 16"
5"	7"	7"	3"
500 lb.	300 lb.	500 lb.	250 lb.
Meshtex	Meshtex	Meshtex	Spacetex /Onyx Vinyl
			n/a
Fits up to 1" wider	Fits 1" narrower up to 1" wider	Fits up to 1" wider	Fits 1" narrower up to 1" wider
SEB1M0010	SEB1M0015	SEB1M0020	SE0001460 SE0001461 SE0001462 SE0001463
INVACARE	INVACARE	INVACARE	MOTION CONCEPTS
MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS	MOTION CONCEPTS

** Slip Covers available for specific back sizes. Additional surcharges apply. Please contact Motion Concepts for details.
†(21", 23", 25"-30") back widths are available as special order only. Please contact Motion Concepts for details.

For more information on Invacare Matrix
Seating Series Products visit
www.invacare.ca

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts

www.motionconcepts.com
info@motionconcepts.com

USA

700 Ensminger Rd., Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada

84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905)-695-0134
Fax: 905-695-0138

© 2018

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
©. All trademarks are
owned by or licensed to
Invacare Corporation unless
otherwise noted. Matrix
is a registered trademark
of Motion Concepts.
Specifications are subject
to change without
notification.
Form No. 17-355C
Rev. 09/18
160955

