

INVACARE®
MyOn™ Junior
Wheelchair



Yes, you can.®

The **Invacare® Myon™ Junior wheelchair** is a lightweight, foldable pediatric wheelchair designed for children aged between 3 and 15 years. MyOn Junior has been developed to match the individual needs of the child and can grow as they grow. Because a child's needs will change the MyOn Junior offers a large range of options to accompany the child in their development and to provide the necessary clinical support. All models of MyOn Junior come in attractive, modern styles, designed to appeal to kids

Highlights

- Width of the seat and backrest can be adjusted, as can the height of the arms and footrests.
- Extensive ranges of size combinations, offering the right fit for all children, from toddlers to teenagers.
- Very suitable for transportation and storage, due to its lightweight design, folding backrest, quick release wheels and detachable arm and legrests.
- The optional stability kit, which consists of rigged out front castors, helps prevent the child from tipping forward when climbing up or down on the footplates.



Specifications

Overall		
Width Open	7.5" plus seat width	
Seat Widths	8"-11"(child)/ 12"-15" (young)	
Seat Depths	8"-10", 10"-12", 12"-14" (child)/ 12"-14", 11"-13", 13-15" (young)	
Front Seat to Floor Heights		
Rear Wheel Size	Seat Slope 0°*, Front STF	
20"	15", 16", 17"	
22"	15", 16", 17", 18"	
*Various combinations can be done. Seat Slope 3° and 5° also available. See charts on current price list for complete details.		
Kneel to Heel Range		
80° Straight Short Swing In/Out	7.5"-13.00"	
80° Straight Long Swing In/Out	12.00"-18.50"	
90° Straight Lift Off	6.75"-12.50"	
Elevating Legrest Lift Off	12.25"-16.25"	
Arm Heights		
T Arm Adj. Height Desk/Junior	6" - 11.25"	
Back Canes		
	Back Canes Degree	Back Canes Heights
Fixed Canes	0°	15"/16"/17"
Adj. Height Canes	7°	14"-16", 15" - 17"
Adj. Angle Canes	7°	14"-18"
Folding Back Canes		15"/16"/17", 16"
Weight	approx. 27 lbs (child), approx. 33lbs (young)	
Weight Capacity	175lbs	

Invacare Canada LP
www.invacare.ca

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2010

All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject to change without notification.

Invacare Canada
Form No.10-789C

