

FREEDOM[®]NXT[™]

LIGHTWEIGHT FOLDING TILT-IN-SPACE MOBILITY



The Freedom NXT is the future of lightweight, folding tilt-in-space mobility. Users and caregivers no longer have to sacrifice compactness for durability or configurability. Conveniently designed to be compact and adjustable, the NXT can grow with you and travel to wherever you need to be.

place your order

1.800.668.5324
phone

1.800.668.5478
fax



FREEDOM NXT

LIGHTWEIGHT FOLDING TILT-IN-SPACE MOBILITY



Lightweight. Strong. Dependable. Compact.

The NXT has proven that strength, durability and reliability does not have to come at the cost of weight, adjustability and portability. It's easy to fold, lift and stow – and in its lightest configuration weights in at only 22 lbs. The robust design has a rigid like feel, and adjusts 2 inches in width and 4 inches in depth to grow with the user while still being able to fit it into the trunk of a small vehicle.

FREEDOM NXT

FEATURES & OPTIONS

Features

The NXT is unique in that it is the first product ever to have a true “specialized” positioning seat frame. With the LEAF & MAPS options, it can fit a seated posture at the hips and at the knees and address foot and ankle positioning challenges.

Even with these special features, there is no impact on the overall size of the NXT.

The NXT is truly the first fully functional positioning frame that remains lightweight, compact and portable and wherever you take it, you can be rest assured it won't let you down!



1. Optional 2" wide casters
2. Optional LEAF frame
3. Optional transport ready frame

Options

- Streamlined and aesthetically pleasing
- Very light 22 lb. transfer weight
- Compact to fit into small spaces
- Folding frame has “rigid-like” feel
- Short, medium & long frame choices
- Rider weight capacity of 225 lbs.
- Easy to adjust for growth
- 21 colours to choose from
- Eleven sizes in 1" increments
- Aluminum front casters
- Universal axle plate
- Patented dual locking safety axle
- Most sizes grow 2" in frame width
- All sizes grow 4" in seat depth
- Back angle adjusts from 80° to 120
- 2", 3" and 4" seat frame extensions
- Anti-rotation seat frame
- Anti-rotation cross brace
- Flared front base frame
- Seat tube leveling guides
- Lightweight front fork
- Foot activated tilt
- Foot activated wheel lock
- Pop Up or Swing Away footrest hangers
- Lower Extremity Accommodation Frame
- Seat Rail Only seat frame option
- Asymmetric seat frame adjustment
- 5° anterior tilt frame setting
- Seating Dynamics Rocker Back option
- Single or dual hand lever activated tilt
- Attendant control hand lever activated tilt, brake and wheel lock
- Length adjustable footrest, extension tubes
- Combination anti-tipper & curb assist leverage bar
- Oxygen tank holder

FREEDOM NXT SPECIALIZED FEATURES & OPTIONS

LEAF in Abduction



LEAF in Windswept Right



LEAF in Windswept Left



MAPS in Plantar Flexion



MAPS in Eversion



MAPS in Inversion



LEAF (Lower Extremity Accomodation Frame):

This patent pending design is a creative solution to complex seating needs. Unique to Freedom Designs, it was designed so that the front frame angle can adjust to accommodate abduction, adduction and windswept leg positioning. Optional on left or right or both, it is paired most ideally with 75 degree or 90 degree pop up swingaway hangers. Easily adjustable with just one allen key.

MAPS (Multi Axis Positioning System):

MAPS (Multi Axis Positioning System): The MAPS footrest was designed to be fully adjustable and help in addressing multiple foot and ankle positioning challenges and alleviates the need for additional positioning accessories. Just one footplate hardware adjustment allows for combinations of dorsiflexion or plantarflexion and inversion or eversion, as well as rotation.

Model	Overall Width (2" increments)	Overall Height	Seat Depth	Seat Width (1" increments)	Seat-to-Floor Height (1" increments)	Tilt Range	Adjustable Back Angle	Tubing Diameter
NXT	+ 8.5", +9.5"	32.5" Min or 44" Max	Short Frame: 12"-16" (18" with extensions) Medium Frame: 14"-18" (21" with extensions) Long Frame: 16"-22" (26" with extensions)	10" to 20"	12.5" - 18" in 1/4" increments	-5° Anterior 45° Posterior	80° to 120°	1"

Back Height	Caster Height	Wheel Height	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty	Total Weight**
20", 23", 26"	3", 4", 5", 6", 7"	12", 16", 18", 20", 22", 24"	225 lbs.	55 lbs. - 65 lbs.	29" x 14" x 39"	5 Years Frame/Crossbases	47 lbs.

Invacare Canada LP
www.invacare.ca

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324
16773, boulevard Hymus
Kirkland, Quebec
H9H 3L4
(800) 363-0137

©2017
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted.
Specifications are subject
to change without notification.
Form No. 17-361C

Rev. 6/18

