

 **TOP END**[®]



Invacare[®] Top End[®] Pro[™] Tennis Wheelchair

For new players, adjustability is key. The **Invacare[®] Top End[®] Pro[™] Tennis Wheelchair** has just that, featuring adjustable front and rear seat heights, center-of-gravity, footrest height, backrest angle/height and anti-tip. Plus, the Pro Tennis Wheelchair is available in two frame height choices: short or tall. The short frame is ideal for children and the tall frame is for adults, which makes this chair a great choice for programs that have multiple users. Adjusts in minutes, super easy to use, lightweight and durable.

Features

- Adjustable front and rear seat heights
- Anti-tip/5th wheel: adjustable height
- Fixed or folding back type
- Adjustable back height
- Adjustable tension back upholstery
- Adjustable center-of-gravity
- Fixed camber 15°, 18° or 20°
- High performance Sun® wheels: 24", 25", 26"
- Anodized handrims
- High pressure clincher tires
- Smart fork caster system
- 72mm front casters and 59mm anti-tip
- Height adjustable platform footrest
- Powder coat paint
- 250 lb. weight capacity

Model #: SPTN



1. Front seat height and center-of-gravity adjustments
2. Rear seat height adjustment
3. Quick release lever
4. Single swivel fixed fifth wheel
5. Precision fork caster system
6. Adjustable footrest system



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(800) 532-8677
Fax: (727) 522-1007

©2018

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols™ and®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. Spinerig
is a registered trademark
of Spinerig Inc. Sun and
Equalizer are registered
trademarks of Sun
Components. Kenda is a
registered trademark of
Kenda Rubber Industrial
Co., Ltd. Specifications are
subject to change without
notification.
Form No. 18-322
181124

Options

- Clamp-on fixed aluminum side guards
- Spinerig® SLX wheels
- 12 Tab handrims
- Kenda® tires
- Click straps
- Frame protectors
- Alignment gauge
- Telescoping anti-tip bar
- Cushion
- Folding Back
- Fold-in side guards