


INVACARE® MATRX®

ELAN®

HEADREST



The **Invacare® Matrx® Elan® Headrest** — Takes head support to the next level. A new elegance in head support and comfort, the Elan Headrest is compact, lightweight and highly adjustable, and compliments a full range of pediatric through adult applications.

Call to order

1.800.668.5324
phone

1.800.668.5478
fax



matrx®
SEATING SERIES

Elan Headrest Hardware Features and Functional Benefits


Elegant – Clean, compact, streamlined aesthetics.

Adjustable – Six points of adjustment – infinite positioning opportunities to meet user needs.

Easy – Single tool set up. Optional tool-less adjustment levers provided for ease of position changes.

Reliable – Multi-disk locking provides superior holding strength.

Sizes – Standard adult (15" post) or mini (11" post) available.



Elan Headrest Pads Features and Functional Benefits

Comfortable – Dual density foam provides immersion, prevents bottoming out.

Durable – Co-molded, dual firmness, HR foam surrounds backing plates, can be easily removed for modification, and are made with Ultra-fresh antimicrobial additive.

Lightweight – 45% Lighter than competitive pads.

Versatile – Standard, occipital and multi-adjustable 4-point support models available.
– Aluminum plates can be formed to desired shape.

Covers – Zippered, moisture-resistant, breathable Startex fabric
– Standard or Infection control cover (reverse Startex) available.



Model Description	Product	Dimensions Width X Height	Startex Cover Model #	Reversed Startex Cover (Infection Control) Model #	Limited Warranty
ESP-6	6" Standard Pad	6" x 3.5"	ESP6	ESP6-IC	2 Years
ESP-10	10" Standard Pad	10" x 5"	ESP10	ESP10-IC	2 Years
ESP-14	14" Standard Pad	14" x 5"	ESP14	ESP14-IC	2 Years
EOP-9	9" Occipital Pad	9" x 7"	EOP9	EOP9-IC	2 Years
EOP-12	12" Occipital Pad	12" x 8"	EOP12	EOP12-IC	2 Years
E4POINT	4 Point Pad	11" x 10"	E4POINT	E4POINT-IC	2 Years
MEHW	Elan Standard Hardware	15" Vertical Post	N/A	N/A	2 Years
MEMH	Elan Mini Hardware	11" Vertical Post	N/A	N/A	2 Years

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com
info@motionconcepts.com

USA
700 Ensminger Rd., Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: 905-695-0138

© 2017
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
©. All trademarks are
owned by or licensed to
Invacare Corporation unless
otherwise noted. Matrix
is a registered trademark
of Motion Concepts.
Specifications are subject
to change without
notification.
Form No. 17-358C

