



# GREAT HYGIENE & MOBILITY SOLUTIONS



Yes, you can.®

# The Great Big Wheel Transport Chair Aluminum



Model no. B10218LR



Model no. B10218LB

## The Great Big Wheel Transport Chair Aluminum

Model no.	B10218LB - 18" Seat Width	Colour: Electric Blue
Model no.	B10220LB - 20" Seat Width	Colour: Electric Blue
Model no.	B10218LR - 18" Seat Width	Colour: Electric Red
Model no.	B10220LR - 20" Seat Width	Colour: Electric Red

The Great Big Wheel Aluminum transport chair delivers a lightweight, compact mobility solution with the added convenience of large 12" rear wheels and dual attendant hand brakes. The Great Big Wheel Aluminum transport chair is personalized for you in two seat width's 18" & 20" with your choice of attractive frame finish Electric Blue and Electric Red. The durable padded nylon seating and padded full length vinyl arm supports ensure comfortable travel wherever life takes you.

### Features

- 2" x 8" and 2" x 12" durable composite spoked wheels with non-marking urethane tires
- Dual braking/wheel lock controls on the handles allow the rear wheels to be slowed or locked in place
- The braking functions and wheel locks operate independently offering maximum control. Standard left and right rear-mounted curb aids, in conjunction with the large 12" wheels, allow the chair to easily overcome curbs and thresholds
- Removable footrests with flip-up footplates offer adjustable comfort and ease of transfer for the user
- Calf Support padded nylon with hook and loop closure
- Easy transport and storage – convenient fold-down back support and folding frame with carry handles allows for simple and easily managed storage and travel
- A rear storage pouch with hook and loop enclosure helps keep personal belongings secure and accessible

### Specifications

Seat Width	18" or 20"
Seat Depth	16"
Back Height	18"
Seat-to-floor height	19"
Footplate-to-seat	15" - 19.5"
Product Weight	19.5 lbs.
User Weight Capacity	250 lbs.
Warranty	1 Year

# The Great Transport Chair Aluminum



Model no. B30118LR



Model no. B30118LB

## The Great Transport Chair Aluminum

Model no.	B30118LB - 18" Seat Width	Colour: Electric Blue
Model no.	B30120LB - 20" Seat Width	Colour: Electric Blue
Model no.	B30118LR - 18" Seat Width	Colour: Electric Red
Model no.	B30120LR - 20" Seat Width	Colour: Electric Red

The Great Transport Chair Aluminium delivers a lightweight, compact mobility solution, packed with standard comfort and convenience features. The Great Transport Chair Aluminium is personalized for you in two seat width's 18" & 20" with your choice of attractive frame finish Electric Blue and Electric Red. The durable padded nylon seating and padded full length vinyl arm supports ensure comfortable travel wherever life takes you.

### Features

- 4" x 8" durable composite spoked wheels with non-marking urethane tires
- Rear wheel locks keep the chair stable for parking and transfer. Standard left and right rear mounted curb aids allow the attendant to easily overcome curbs and thresholds
- Removable footrests with flip-up footplates offer adjustable comfort and ease of transfer for the user
- Calf Support padded nylon with hook and loop closure
- Easy transport and storage – convenient fold-down back support and folding frame with carry handles allows for simple and easily managed storage and travel
- A rear storage pouch with hook and loop enclosure helps keep personal belongings secure and accessible

### Specifications

Seat Width	18" or 20"
Seat Depth	16"
Back Height	18"
Seat-to-floor height	19"
Footplate-to-seat	15" - 19.5"
Product Weight	19.5 lbs.
User Weight Capacity	250 lbs.
Warranty	1 Year

# The Great Big Wheel Transport Chair Steel



Model no. B10218

## The Great Big Wheel Transport Chair Steel

Model no. B10218 - 18" Seat Width      Colour: Silver Vein  
Model no. B10220 - 20" Seat Width      Colour: Silver Vein

The Great Big Wheel Transport Chair Steel delivers a compact mobility solution with durable steel construction and 300 lbs. weight capacity. In addition the added convenience of large 12" rear wheels and dual attendant hand brakes. The Great Big Wheel Transport Chair Steel is personalized for you in two seat width's 18" & 20" with an attractive resilient silver vein finish. The durable padded nylon seating and padded full length vinyl arm supports ensure comfortable travel wherever life takes you.

### Features

- 2" x 8" and 2" x 12" durable composite spoked wheels with non-marking urethane tires
- Dual braking/wheel lock controls on the handles allow the rear wheels to be slowed or locked in place
- The braking functions and wheel locks operate independently offering maximum control. Standard left and right rear-mounted curb aids, in conjunction with the large 12" wheels, allow the chair to easily overcome curbs and thresholds
- Removable footrests with flip-up footplates offer adjustable comfort and ease of transfer for the user
- Calf Support padded nylon with hook and loop closure
- Easy transport and storage – convenient fold-down back support and folding frame with carry handles allows for simple and easily managed storage and travel
- A rear storage pouch with hook and loop enclosure helps keep personal belongings secure and accessible

### Specifications

Seat Width	18" or 20"
Seat Depth	16"
Back Height	18"
Seat-to-floor height	19"
Footplate-to-seat	15" - 19.5"
Product Weight	24.5 lbs.
User Weight Capacity	300 lbs.
Warranty	1 Year

# The Great Transport Chair Steel



Model no. B30118

## The Great Transport Chair Steel

Model no. B30118 - 18" Seat Width      Colour: Silver Vein  
Model no. B30120 - 20" Seat Width      Colour: Silver Vein

The Great Transport Chair Steel delivers a compact mobility solution with durable steel construction and 300lb user weight capacity. The Great Transport Chair Steel is personalised for you in two seat width's 18" & 20" with an attractive resilient silver vein finish. The durable padded nylon seating and padded full length vinyl arm supports ensure comfortable travel wherever life takes you.

### Features

- 4" x 8" durable composite spoked wheels with non-marking urethane tires
- Rear wheel locks keep the chair stable for parking and transfer. Standard left and right rear mounted curb aids allow the attendant to easily overcome curbs and thresholds
- Removable footrests with flip-up footplates offer adjustable comfort and ease of transfer for the user
- Calf Support padded nylon with hook and loop closure
- Easy transport and storage – convenient fold-down back support and folding frame with carry handles allows for simple and easily managed storage and travel
- A rear storage pouch with hook and loop enclosure helps keep personal belongings secure and accessible

### Specifications

Seat Width	18" or 20"
Seat Depth	16"
Back Height	18"
Seat-to-floor height	19"
Footplate-to-seat	15" - 19.5"
Product Weight	22.5 lbs.
User Weight Capacity	300 lbs.
Warranty	1 Year

# The Great Shower Chair With and Without Back



Model no. B231A



Model no. B235A

## The Great Shower Chair

Model no. B231A - The Great Shower Chair

Model no. B235A - The Great Shower Chair With Back

The Great Shower Chairs with and without back are different and better because it was designed for you. Simple to assemble and quick to take apart. Shower, travel, storage, you choose when and where without the need for tools and fasteners that strip.

The Great Shower chair B235A with back is finished with 4x 2.5" suction tips with composite reinforcement inserts providing an additional level of stability, durability and corrosion resistance.

## Features

- Shower in comfort enjoying the large contoured seating area
- The one-piece injection molded polycarbonate seat supports up to 500lbs and is stain resistant
- The adjustable aluminum supporting frame ensures perfect seat height and will never rust

## Specifications

	B231A	B235A
Assembly	Tool-free/Push button	Tool-free/Push button
Seat Height Adjustment	14.5" - 18.5"	14.5" - 18.5"
Seat Dimension	17" W x 14.5" D	17" W x 14.5" D
Back Dimension	N/A	14.5" W x 9.5" H
Seat To Top Of Back	N/A	12"
Outside Width At Base	18.5"	18.5"
Outside Depth At Base	16"	16"
Product Weight	4.8 lbs.	6.1 lbs.
User Weight Capacity	500 lbs.	500 lbs.
Warranty	13 Months	13 Months

# The Great Shower Stool



Model no. B145A

## The Great Shower Stool

Model no. B145A - The Great Shower Stool

The Great Shower Stool. No need to compromise on quality and features. The Great Shower Stool contains a durable aluminum frame, supporting up to 300 lbs. and ample seating area. It is equipped with integrated carrying and positioning handles.

### Features

- Tool-free assembly ensures set-up and also makes disassembly a breeze for storage or travel
- The push button adjustable legs provide the perfect seat height
- The Great Shower Stool is rust free, finished with 4" x 2.5" suction tips with composite reinforcement inserts, providing an additional level of stability, durability and corrosion resistance

### Specifications

Assembly	Tool-free/Push button
Seat Height Adjustment	13" - 19"
Seat Dimension	12.5" Diameter
Outside Width At Base	19"
Product Weight	3.0 lbs.
User Weight Capacity	300 lbs.
Warranty	13 Months

# The Great Transfer Bench



Model no. B281C



Model no. B281CUP

## The Great Transfer Bench

Model no. B281C – The Great Tool-Free Transfer Bench

Model no. B281CUP – Great Tool Free Transfer Bench Large Suction Legs Accessory, Pair

The Great Tool-Free Transfer Bench makes transfers in and out of your bath tub comfortable and secure. With a quick and simple slide transfer, a user eliminates the potentially dangerous step over the tub wall.

Lightweight and portable, the blue Great tool free transfer bench is the perfect hygiene companion at home or away.

The Great Tool Free Transfer large suction tips feature 2.5" suction tips for added stability. The B281CUP are an add-on option for the B281C Great Tool Free Transfer Bench. These legs can be used inside, outside or on both sides of the tub to provide maximum stability. The B281C tips can be released by pulling up on the marked release tabs on the cups.

## Features

- The comfortable seating area provides a back support and handle to aid in secure positioning and safety
- The durable aluminum frame is corrosion proof and height adjustable, ensuring the perfect seat height to accommodate most standard tubs
- The B281C is finished with four latex-free suction tips providing stability and protection to tub finishes

## Specifications

Assembly	Tool-free/Push button
Seat Height Adjustment	16.9" – 20.7"
Seat Dimension	26.4" W x 16" D
Outside Width At Base	28.9"
Outside Depth At Base	16.7"
Product Weight	8.4 lbs.
User Weight Capacity	250 lbs.
Warranty	13 Months



# The Great Round and Long Seat Risers



Model no. B255A



Model no. B251A

## The Great Round and Long Seat Risers

Model no. B255A - The Great Round Seat Riser  
Model no. B251A - The Great Long Seat Riser (Elongated)

All the comfort, support and stability of traditional risers without clumsy attachments and undesirable appearance. Feel right at home with the extra support you need. The great risers mount securely to any toilet maintaining the use of the existing seat and lid, without reducing seat opening. Cleaning is a breeze thanks to the patented interlock mounting system. Interlock enables simple and easy removal for cleaning anytime.

### Features

- Simple installation
- Two piece design enables easy cleaning and infection control
- Safe and secure mounting (hardware included)
- Clean and simple aesthetics
- Fits residential and commercial round toilets
- Riser height - 3" plus height of existing seat and lid total rise 4" - 4.5"

### Specifications

	B255A: Round	B251A: Elongated
Inside Opening	9.5" x 11.5"	19.5" x 14.0"
Outside Dimension	14.5" x 17.75"	14.25" x 20.25"
Product Weight	2.70 lbs.	2.75 lbs.
User Weight Capacity	300 lbs.	300 lbs.
Warranty	13 Months	13 Months

# The Great Toilet Safety Frame



Model no. B1392

## The Great Toilet Safety Frame

---

Model no. B1392 - The Great Toilet Safety Frame

The Great Toilet Safety Frame ensures stability and comfort for transfers to and from the toilet. The B1392 mounts securely and easily using the existing seat hardware and fits all residential and commercial toilets. The Great Toilet Safety Frame provides superior level of stability using the toilet and the floor to distribute weight ensuring a safe transfer for users.

Cleaning and care are simple, and the sturdy aluminum frame will never rust. The B1392 can be used in conjunction with the Invacare 1300, 1301 and 1302 raised toilet seats.

## Features

---

- Simple installation
- Safe and secure mounting (with existing seat hardware)
- Clean and simple aesthetics
- Fits residential and commercial round toilets
- The aluminum and composite construction are easy to clean and will never rust

## Specifications

---

Width Inside Frame	18.5" - 20.5"
Width Outside Frame	20.5 - 22.5"
Height At Handles	25" - 28"
Product Weight	5.5 lbs.
User Weight Capacity	300 lbs.
Warranty	13 Months

---

# The Great Bed Support Rail



Model no. B1700

## The Great Bed Support Rail

Model no. B1700 - The Great Bed Support Rail

The Great Bed Support Rail will fit onto almost any traditional mattress and box spring bed.

The B1700 is installed tool-free, positioned underneath the mattress and secured using the provided strap. The adjustable height handle allows it to work with mattress' of varying sizes and accommodating the perfect height for you. The Great Bed Support Rail is designed with a chrome finish and foam handle.

The unique rigid construction when assembled, eliminates flex and provides superior strength not found in similar commercially available products.

## Features

- Tool-free and simple installation
- Fits traditional mattress and box spring bed
- The chrome finish and foam handle adapt with any decor

## Specifications

Width	20"
Length	35.5"
Height	18" - 22"
Product Weight	3.5 lbs.
User Weight Capacity	265 lbs.
Warranty	13 Months



*Yes, you can.®*

**Invacare Canada LP**

[www.invacare.ca](http://www.invacare.ca)

570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

All rights reserved.  
Trademarks are identified  
by the symbols <sup>TM</sup>, <sup>SM</sup> and <sup>®</sup>.  
All trademarks are owned  
by or licensed to Invacare  
Corporation unless  
otherwise noted.

Specifications are subject to  
change without notification.

© 2019  
Invacare Canada  
Form No. 17-351C  
Rev 03/19