



INVACARE® TOP END®  
**FORCE™ K**  
HANDCYCLE



The **Invacare® Top End® Force™ K Handcycle** is a kneeling position handcycle available in two custom, made-to-measure aerodynamic designs. The Invacare Top End Force K Kneeler is used by the most able athletes, while the Top End Force K Double Amputee Kneeler is custom designed for the double above-the-knee amputee. What makes the Force K Handcycle different from the rest of the Invacare Top End Force Series handcycles? The Force K Handcycle is a kneeling-only, non-adjustable position handcycle that requires dialed-in measurements and engineering drawings—all of which are easy to find on the website: [www.topendwheelchair.com](http://www.topendwheelchair.com). Can you force the pace? Get the Force K Handcycle.

## Features

7005 aluminum custom welded frame  
7005 aluminum aerodynamic fork  
Breathable mesh upholstery  
Seat and back cushions  
Tapered crank arms  
Frictionless pedal bearings  
Aluminum ergonomic hand pedals  
SRAM® 30 speed drive train  
Chainrings: 53-39 or 48-39-24  
Cassettes: SRAM® 11-34

Tektro® Caliper Brake  
Fork mounted locking second brake  
Steering dampener  
26" high performance wheels  
3°, 6°, 9° or 11° of camber  
Threaded axles  
Drafting bumper  
Parking brake  
Safety flag  
Chain guard

Many paint choices  
Reflective stickers  
250 lb. weight capacity

## Model #: FRCK



1. Top End Carbon Wheel Upgrade is stiff, light and aerodynamic
2. SRAM components allow shifting and braking at your fingertips
3. Adjustable depth seat upholstery in short position
4. Adjustable depth seat upholstery in long position

**Invacare Corporation**  
[www.invacare.com](http://www.invacare.com)

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

### Top End

[www.topendwheelchair.com](http://www.topendwheelchair.com)  
Sports and Recreation Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(727) 522-8677  
(800) 532-8677  
Fax: (727) 522-1007

©2016

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and ®. All trademarks are  
owned by or licensed to  
Invacare Corporation unless  
otherwise noted. SRM is a  
trademark of SRM Service  
Center Inc. USA. SRAM is a  
registered trademark of SRAM  
Corp. Tektro is a registered  
trademark of Tektro USA.  
Specifications are subject to  
change without notification.  
Form No. 12-451  
Rev. 03/16  
151184

## Options

Carbon fiber wheels  
Component, brake and cassette upgrades  
SRM™ Power Meter measuring tool integrated  
directly into Invacare Top End cranks

See the Force K Handcycle order form for all of the  
great options available!

## Frame Styles:

The Force K Kneeler is designed for the most able  
athletes with a custom knee/seat bucket.

The Force K Double Amputee Kneeler is designed  
for the double above-the-knee amputee utilizing a  
custom cage.

