

 **TOP END**®



Invacare® Top End® Pro™ Basketball Wheelchair

The **Invacare® Top End® Pro™ Basketball Wheelchair** delivers custom-built performance with adjustments that take only minutes. Perfect for competitors still dialing in their game.

Features

- Custom aluminum frame
- Offensive 4 bend wing
- Adjustable front and rear seat heights
- Anti-tip/fifth wheel: Adjustable height
- Fixed or folding back type
- Adjustable back height
- Adjustable tension back upholstery
- Adjustable Center-of-Gravity
- Rear seat height adjustment range is 6" and the frame is available in a tall or short variety to accommodate every aspiring athlete
- Fixed camber: 15°, 18° or 20°
- High performance Sun® wheels
- Anodized handrims
- High pressure clincher tires
- Smart fork caster system
- 72mm casters / 59mm rear anti-tip
- Height adjustable platform footrest
- Powder coat paint
- 250 lb. weight capacity
- Quick release lever allows the player to change his or her rear seat height in just minutes

Model #: SPBB



1. Front seat height and center-of-gravity adjustments
2. Rear seat height adjustment
3. Quick release lever
4. Single anti-tip
5. Precision fork caster system
6. Adjustable footrest system



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(800) 532-8677
Fax: (727) 522-1007

©2018

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols™ and®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. Spinerigy
is a registered trademark
of Spinerigy Inc. Sun and
Equalizer are registered
trademarks of Sun
Components. Kenda is a
registered trademark of
Kenda Rubber Industrial
Co., Ltd. Specifications are
subject to change without
notification.
Form No. 18-318
181120

Options

- Clamp-on fixed aluminum side guards
- Spinerigy® SLX wheels
- 12 tab handrims
- O-Ring positioning belt
- Close mount/cut pushrim tab
- Click straps
- Frame protectors
- Alignment gauge
- Double swivel anti-tip bar
- Cushion: soft or hard
- Folding back

 **TOP END®**